OsteoBiol[®]
by Tecnos

YEARS
25
ANNIVERSARY

Bone Grafting Materials

REGENERATION SCIENCE

INSPIRED BY NATURE

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
Giuseppe Oliva MD
R&D Director
Tecnoss S.r.l.

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(1) Figueiredo M, Henriques J, Martins G, Guerra F, Judas F, Figueiredo H
Physicochemical characterization of biomaterials commonly used in dentistry as bone substitutes - comparison with human bone
Journal of Biomedical Materials Research Part B: Applied Biomaterials, 2010 Feb; 92(2):409-19

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(2) Nannmark U, Sennerby L

The bone tissue responses to prehydrated and collagenated cortico-cancellous porcine bone grafts: a study in rabbit maxillary defects
Clinical Implant Dentistry and Related Research, 2008 dec;10(4):264-70

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Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy
Clinical Oral Investigations, 2018 Jan;22(1):505-513

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
(5) Brunelli G, Sollazzo V, Carinci F, Palmieri A, Girardi A, Monguzzi R
OsteoBiol® influences osteogenic differentiation of adipose derived stem cells
European Journal of Inflammation, 2011, Vol. 9, no. 3 (S), 103-107

(6) Jeanneau C, Le Fournis C, About I
Xenogeneic bone filling materials modulate mesenchymal stem cell recruitment: role of the complement c5a
Clinical oral investigations; 2019 oct 23


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5

Fig. 6

Fig. 7


Fig. 8

Fig. 9

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
Wide diameter immediate post-extractive implants vs delayed placement of normal-diameter implants in preserved sockets in the molar region: 1-year post-loading outcome of a randomised controlled trial
European Journal of Oral Implantology, 2017;10(3):263-278

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


(8) Barone A, Orlando B, Tonelli P, Covani U
Survival rate for implants placed in the posterior maxilla with and without sinus augmentation: a comparative cohort study
Journal of Periodontology, 2011 Feb; 82(2):219-26

(9) Scarano A, Piattelli A, Assenza B, Quaranta A, Perrotti V, Piattelli M, Iezzi G
Porcine bone used in sinus augmentation procedures: a 5-year retrospective clinical evaluation
Journal of Oral and Maxillofacial Surgery, 2010 Aug;68(8):1869-73

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

BIBLIOGRAFIA

- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (4) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42
- (7) MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR
MATERIALS, 2017 SEP 8;10(9)
- (8) CHECCHI V, FELICE P, ZUCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278
- (9) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- (10) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (11) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRA-BONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (12) ROMBOUTS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

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

BIBLIOGRAFIA

- (1) CORREIA F, GOUVEIA SA, POZZA DH, FELINO AC, FARIA-ALMEIDA R.
A RANDOMIZED CLINICAL TRIAL COMPARING IMPLANTS PLACED IN TWO DIFFERENT BIOMATERIALS USED FOR MAXILLARY SINUS AUGMENTATION.
MATERIALS. 2023; 16(3):1220.
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVESTIG, 2017 2018 JAN;22(1):505-513
- (4) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870
- (5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (6) SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- (7) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (8) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (9) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (10) BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759
- (11) PISTILLI R, CANULLO L, PESCE P, PISTILLI V, CAPONIO C, SBRICOLI L
GUIDED IMPLANT SURGERY AND SINUS LIFT IN SEVERELY RESORBED MAXILLAE: A RETROSPECTIVE CLINICAL STUDY WITH UP TO 10 YEARS OF FOLLOW-UP
J DENT, 2022 APR 21;121:104137. ONLINE AHEAD OF PRINT
- (12) CORREIA F, GOUVEIA SA, POZZA DH, FELINO AC, FARIA-ALMEIDA R.
A RANDOMIZED CLINICAL TRIAL COMPARING IMPLANTS PLACED IN TWO DIFFERENT BIOMATERIALS USED FOR MAXILLARY SINUS AUGMENTATION
MAT. 2023; 16(3):1220

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
BIBLIOGRAFIA

- (1) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005 JUN;54(6):351-62
- (2) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006;2: 99-106
- (3) ROMASCO T, TUMEDEI M, INCHINGOLO F, PIGNATELLI P, MONTESANI L, IEZZI G, PETRINI M, PIATTELLI A, DI PIETRO N
A NARRATIVE REVIEW ON THE EFFECTIVENESS OF BONE REGENERATION PROCEDURES WITH OSTEOBIOL® COLLAGENATED PORCINE GRAFTS: THE TRANSLATIONAL RESEARCH EXPERIENCE OVER 20 YEARS
J FUNCT BIOMATER, 2022 AUG 18;13(3):121
- (4) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (5) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (6) SANTAGATA M, GUARNIELLO L, TARTARO G
A MODIFIED EDENTULOUS EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N.:114-9
- (7) MALTAGLIATI A, ANGIERO F, FERRANTE F, BLASI S, OTTONELLO A
ALVEOLAR RIDGE EXPANSION BY IMPLANTS WITH OSTEODISTRACTION: A CLINICAL REPORT
DENT ORAL CRANIOFAC RES. 2016, VOLUME 2, ISSUE 6
- (8) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

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
SINUS LIFT
mini-rialzo crestale

BIBLIOGRAFIA

- (1) LOMBARDI T, LAMAZZA L, BERNARDELLO F, ZIĘTEK G, STACCHI C, TROIANO G
CLINICAL AND RADIOGRAPHIC OUTCOMES FOLLOWING TRANSCRESTAL MAXILLARY SINUS FLOOR ELEVATION WITH INJECTABLE XENOGENOUS BONE SUBSTITUTE IN GEL FORM: A PROSPECTIVE MULTICENTER STUDY
INT J IMPLANT DENT. 2022 JUL 22;8(1):32
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL. 2010;36(6):485-9
- (3) DEVELIOGLU H, OZCAN G, GULTEKIN SE, SENGUVEN B, YILDIRIM A
THE SHORT-TERM EFFECTS OF VARIOUS XENOGRAFTS ON BONE HEALING IN RATS CRANIAL DEFECTS
BIOMEDICAL RESEARCH 2018;29(8); 1598-1602
- (4) COMUZZI L, TUMEDEI M, PIATTELLI A, TARTAGLIA G, DEL FABBRO M
RADIOGRAPHIC ANALYSIS OF GRAFT DIMENSIONAL CHANGES IN TRANSCRESTAL MAXILLARY SINUS AUGMENTATION: A RETROSPECTIVE STUDY
MATERIALS (BASEL). 2022 APR 19;15(9):2964
- (5) BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT. 2008 JUN;28(3):283-9
- (6) STACCHI C, SPINATO S, LOMBARDI T, BERNARDELLO F, BERTOLDI C, ZAFFE D, NEVINS M
MINIMALLY INVASIVE MANAGEMENT OF IMPLANT-SUPPORTED REHABILITATION IN THE POSTERIOR MAXILLA, PART II. SURGICAL TECHNIQUES AND DECISION TREE
INT J PERIODONTICS REST DENT. 2020 MAY/JUN;40(3):E95-E102
- (7) COMUZZI L, TUMEDEI M, PETRINI M, ROMASCO T, LORUSSO F, DE ANGELIS F, PIATTELLI A, TATULLO M, DI PIETRO N
CLINICAL AND RADIOLOGICAL EVALUATION OF A SELF-CONDENSING BONE IMPLANT IN ONE-STAGE SINUS AUGMENTATION: A 3-YEAR FOLLOW-UP RETROSPECTIVE STUDY
INT. J. ENVIRON. RES. PUBLIC HEALTH 2023, 20, 2583
- (8) COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS. 2008 SEP-OCT;23(5):841-6
- (9) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A HEMINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT. 2009 FEB;29(1):59-67
- (10) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES. 2010 JUN 1;12(2):161-3
- (11) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES. 2012 OCT;14(5):746-58

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
Fonte: Tecnos® Dental Media Library

BIBLIOGRAFIA

- (1) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- (2)ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (3) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107
- (4) CAKIR M, KARACA IR, AYSEGÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION
CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85
- (5) KOLMAS J, SZWAJA M, KOLODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENIC BIOMATERIALS USED AS BONE SUBSTITUTES
J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41
- (6) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U
CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL
J CLIN PERIODONTOL, 2017 FEB;44(2):204-214
- (7) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- (8) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- (9) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- (10) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76
- (11) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (12) MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A
CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL
CLIN IMPL DENT RELAT RES, 2018 Dec;20(6):906-914

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
Fonte: Tecnos® Dental Media Library

BIBLIOGRAFIA

- ((1) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870
- ((2) ALQUTUB M, MUKHTAR A, ALALI Y, VOHRA F, ABDULJABBAR T
OSTEOGENIC DIFFERENTIATION OF PERIODONTAL LIGAMENT STEM CELLS SEEDED ON EQUINE-DERIVED XENOGRAFT IN OSTEOGENIC GROWTH MEDIA
MEDICINA (KAUNAS) 2022 OCT 25;58(11):1518
- ((3) ESPOSITO M, BARAUSSE C, PISTILLI R, PIATELLI M, DI SIMONE S, IPPOLITO DR, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54
- ((4) FELICE P, BARAUSSE C, BARONE A, ZUCCHELLI G, PIATELLI M, PISTILLI R, IPPOLITO DR, SIMION M
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480
- ((5) BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447
- ((6) MARCONCINI S, COVANI U, GIAMMARINARO E, VELASCO-ORTEGA E, DE SANTIS D, ALFONSI F, BARONE A
CLINICAL SUCCESS OF DENTAL IMPLANTS PLACED IN POSTERIOR MANDIBLE AUGMENTED WITH INTERPOSITIONAL BLOCK GRAFT: 3-YEAR RESULTS FROM A PROSPECTIVE COHORT CLINICAL STUDY
J ORAL MAXILLOFAC SURG. 2019 FEB;77(2):289-298
- ((7) BERNARDI S, GATTO R, SEVERINO M, BOTTICELLI G, CARUSO S, RASTELLI C, LUPI E, ROIAS AQ, IACOMINO E, FALISI G
SHORT VERSUS LONGER IMPLANTS IN MANDIBULAR ALVEOLAR RIDGE AUGMENTED USING OSTEOGENIC DISTRACTION: ONE-YEAR FOLLOW-UP OF A RANDOMIZED SPLIT-MOUTH TRIAL
J ORAL IMPLANTOL, 2018 JUN;44(3):184-191
- ((8) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR-APR;33(2):159-66
- ((9) FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
IMPLANT DENT, 2012 JUN;21(3):190-5

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
BIBLIOGRAFIA

- (1) CABALLÉ-SERRANO J, MUNAR-FRAUA A, DELGADO L, PÉREZ R, HERNÁNDEZ-ALFARO F
PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER MEMBRANES FOR BONE REGENERATION
J MECH BEHAV BIOMED, 2019;97:13-20
- (2) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (3) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58
- (4) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT AND RELAT RES, 2013 OCT;15(5):707-13
- (5) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (6) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85
- (7) POLIS-YANES C, CADENAS-SEBASTIÁN C, GUAL-VAQUÉS P, AYUSO-MONTERO R, MARI-ROIG A, LÓPEZ-LÓPEZ J
GUIDED BONE REGENERATION OF AN ATROPHIC MAXILLA USING HETEROLOGOUS CORTICAL LAMINA
CASE REP DENT., 2019; 5216362
- (8) SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F
A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592
- (9) SCHUH PL, WACHTEL H, BEUER F, GOKER F, DEL FABBRO M, FRANCIETTI L, TESTORI T
MULTI-LAYER TECHNIQUE (MLT) WITH PORCINE COLLAGENATED CORTICAL BONE LAMINA FOR BONE REGENERATION PROCEDURES AND IMMEDIATE POST-EXTRACTION IMPLANTATION IN THE ESTHETIC AREA: A RETROSPECTIVE CASE SERIES WITH A MEAN FOLLOW-UP OF 5 YEARS
MATERIALS, 2021 SEP 9;14(18):5180
- (10) GHALLAB NA, ELASKARY A, ELSABAGH H, TOUKHY AE, ABDELRAHMAN H, EL-KIMARY G
A NOVEL ATRAUMATIC EXTRACTION TECHNIQUE USING VESTIBULAR SOCKET THERAPY FOR IMMEDIATE IMPLANT PLACEMENT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
ORAL MAXILLOFAC SURG, 2022 JUN 20. ONLINE AHEAD OF PRINT
- (11) LOPEZ MA, PASSARELLI PC, NETTI A, D'ADDONA A, CARINCI F, WYCHOWANSKI P, CECCHETTI F
PORCINE CORTICAL BONE LAMINA AS A PREDICTABLE TECHNIQUE FOR GUIDED BONE REGENERATION: HISTOMORPHOMETRIC AND RADIOGRAPHIC EVALUATION
APPLIED SCIENCES 2022, 12(20), 10285

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





- **A)** Barone A et al. - **Flap versus flapless procedure for ridge preservation in alveolar extraction sockets: a histological evaluation in a randomized clinical trial**
Clinical Oral Implants Research, 2015 Jul;26(7):806-13
- **B)** Giuliani A et al. - **Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy**
Clinical Oral Investigations, 2018 Jan;22(1):505-513
- **C)** Barone A et al. - **Clinical and histological changes after ridge preservation with two xenografts: preliminary results from a multicenter randomized controlled clinical trial**
Journal of Clinical Periodontology, 2017 Feb;44(2):204-214
- **D)** Crespi R et al. - **Corticancellous porcine bone in the healing of human extraction sockets: combining histomorphometry with osteoblast gene expression profiles in vivo**
Int Journal of Oral and Maxillofacial Implants, 2011 Jul - Aug; 26(4):866-72
- **E)** Crespi R et al. - **Comparison of magnesium-enriched hydroxyapatite and porcine bone in human extraction socket healing: a histologic and histomorphometric evaluation**
Int Journal of Oral and Maxillofacial Implants, 2011 Sep-Oct;26(5):1057-62
- **F)** Barone A et al. - **Xenograft versus extraction alone for ridge preservation after tooth removal: a clinical and histomorphometric study**
Journal of Periodontology, 2008 Aug; 79(8):1370-7

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- **A)** Cassetta M et al. - **Bone formation in sinus augmentation procedures using autologous bone, porcine bone, and a 50 : 50 mixture: a human clinical and histological evaluation at 2 months**
Clinical Oral Implants Research, 2015 Oct; 26(10):1180-4
- **B)** Scarano A et al. - **Maxillary sinus augmentation in humans using cortical porcine bone: a histological and histomorphometrical evaluation after 4 and 6 months**
Clinical Implant Dentistry and Related Research, 2011 Mar; 13(1):13-18
- **C)** Orsini G et al. - **Histologic and ultrastructural analysis of regenerated bone in maxillary sinus augmentation using a porcine bone-derived biomaterial**
Journal of Periodontology, 2006 Dec; 77(12):1984-90
- **D)** Iezzi G et al. - **Comparative histological results of different biomaterials used in sinus augmentation procedures: a human study at 6 months**
Clinical Oral Implants Research, 2012 Dec;23(12):1369-76
- **E)** Tanaka K et al. - **Sinus floor elevation and antrostomy healing: a histomorphometric clinical study in humans**
Implant dentistry, 2019 Dec; 28(6):537-542

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

BIBLIOGRAFIA

- (1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (2) KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
BMC ORAL HEALTH, 2017 OCT 5;17(1):126
- (3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (5) GIULIANI A, IEZZI G, MAZZONI S, PIATELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICRO TOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVEST, 2017 2018 JAN;22(1):505-513
- (6) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6
- (7) SCARANO A, PIATELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (8) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- (9) BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY
CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112
- (10) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRA BONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

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
(3) TALLARICO M, XHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM
SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):263-275

(4) BALDI N, BUTI J, MENSIS M, ALFONSI F, CINQUINI C, TONELLI P, BARONE A
XENOGENIC DERMAL MATRIX VERSUS AUTOLOGOUS CONNECTIVE TISSUE GRAFT VERSUS NO GRAFT AT ABUTMENT CONNECTION FOR IMPROVING AESTHETICS: 6-MONTH OUTCOMES OF A RANDOMISED CONTROLLED TRIAL
CLINICAL TRIALS IN DENTISTRY, 2020;02(2):49-62

(5) FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESSENCE INT, 2014 NOV-DEC;45(10):853-60

(6) FISCHER K R, TESTORI T, WACHTEL H, MÜHLEMANN S, HAPPE A, DEL FABBRO M
SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATEDPORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES
CLIN IMPLANT DENT RELAT RES., 2019;1-8

(7) VERARDI S, ORSINI M, LOMBARDI T, AUSENDA F, TESTORI T, PULICI A, OREGLIA F, VALENTE NA, STACCHI C
COMPARISON BETWEEN TWO DIFFERENT TECHNIQUES FOR PERI-IMPLANT SOFT TISSUE AUGMENTATION: PORCINE DERMAL MATRIX GRAFT VS. TENTING SCREW
J PERIODONTOL. 2020; ACCEPTED, IN PUBLICATION

(8) FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES
QUINTESSENCE INT, 2013;44(3):243-6

(9) MATOH U, PETELIN M, GASPERSIC R
SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESSIONS
INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

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
storia del marchio **OsteoBiol**[®] by Tecnos

INNOVAZIONE

A close-up photograph of a scientist wearing a white lab coat, a white surgical cap, and a white face mask. The scientist is looking through the eyepieces of a white and black microscope. The background is a blurred laboratory setting with various pieces of equipment. The word "INNOVAZIONE" is overlaid in large, bold, blue capital letters on the left side of the image.

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BIBLIOGRAFIA

- (1) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (2) BARONE A, ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS
CURR PHARM BIOTECHNOL, 2017;18(1):64-75
- (3) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (4) BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73
- (7) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (8) CASSETTA M, RICCI L, IEZZI G, DELLAQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (9) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (10) JEANNEAU C, LE FOURNIS C, ABOUT I
XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

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BIBLIOGRAFIA

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

(3) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(4) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6

(5) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

(6) BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HYSTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9

(7) ROMBOULTS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

(8) JEANNEAU C, LE FOURNIS C, ABOUT I
XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

(9) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

(10) MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR MATERIALS, 2017 SEP 8;10(9)

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BIBLIOGRAFIA

- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (3) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (4) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR-APR;33(2):159-66
- (5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017,18,33-44
- (6) MIZUNO M, FUJISAWA R, KUBOKI Y
TYPE I COLLAGEN-INDUCED OSTEOBLASTIC DIFFERENTIATION OF BONE-MARROW CELLS MEDIATED BY COLLAGEN-A2B1 INTEGRIN INTERACTION
J CELL PHYSIOL. 2000 AUG;184(2):207-13
- (7) HSU FY, CHUEH SC, WANG YJ
MICROSPHERES OF HYDROXYAPATITE/RECONSTITUTED COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH
BIOMATERIALS 1999, 20:1931-1936
- (8) ABDELGAWAD ME, SØE K, ANDERSEN TL, MERRILL DM, CHRISTIANSEN P, KJÆRGAARD-ANDERSEN P, DELAISSE JM
DOES COLLAGEN TRIGGER THE RECRUITMENT OF OSTEOBLASTS INTO VACATED BONE RESORPTION LACUNAE DURING BONE REMODELING?
BONE, 2014 OCT;67:181-8
- (9) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907
- (10) JEANNEAU C, LE FOURNIS C, ABOUT I
XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

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

BIBLIOGRAFIA

- (1) GRIFFITH LG, NAUGHTON G
TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES
SCIENCE 2002, 295:1009-14
- (2) REDDI AH
MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS
TISSUE ENG 2000, 6(4):351-59
- (3) NAKASHIMA N, REDDI AH
THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING
NAT BIOTECHNOL 2003, 9:1025-32
- (4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE
ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS
J BIOMED BIOTECHNOL 2004, 1:24-34
- (5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107

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(1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70.

(2) CHECCHI V, FELICE P, ZUCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278

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(1) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN. ORAL INVESTIG., 2017 2018 JAN;22(1):505-513

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
REGENERATION SCIENCE

INSPIRED BY NATURE

and with
SUBSTITUTE
mp2™ Gen-Di
The Best Evidence

LASL

and with
SUBSTITUTE
mp2™ Gen-Di

LA2

and with
SUBSTITUTE
mp2™ Gen-Di
EMBRANE
mp2™ Evolution

ALR

novation:

Lateral Access Sinus Lift
200
Michele Cosentino
Laura Ricci
Giovanni Iuzzi
Sabrina Calassio
Adriano Pignelli
Vittoria Ferraro

OSSEOUS GRAFT
OsteoBiol®
Gen-Di
Tel: +39 02 22174216

Abstract text describing the study on lateral access sinus lift.

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CASSETTA M, CALASSO S, VOZZA I, DELL'AQUILA D
REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY
EUR J IMPLANT PROSTHODONTICS, 2005;3(1):133-144

ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005 JUN;54(6):351-62

BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25

RINNA C, UNGARI C, SALTAREL A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
J CRANIOFAC SURG, 2005 NOV; 16(6):968-72

BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106

BARONE A, SANTINI S, SBORDONE L, CRESPI R, COVANI U
A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5

COVANI U, BARONE A, CORNELINI R, CRESPI R
CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL
J PERIODONTOL, 2006 APR;77(4):722-7

ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90

TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20 (1 SUPPL 1):87-91

BARONE A, COVANI U
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON- VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS
J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46

DEL CORSO M
SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES
COSMETIC DENTISTRY, 2008, 3:16-20

BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L, GHERLONE E, COVANI U
OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS
CLIN ORAL IMPLANTS RES, 2008 MAY;19(5):511-5

BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9

BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
J PERIODONTOL, 2008 AUG;79(8):1370-7

COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6

CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
J OSSEOINTEGR, 2009;1(1):35-40

COVANI U, MARCONCINI S, CRESPI R, BARONE A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
J ORAL IMPLANTOL, 2009; 35(4):189-95

FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010FEB; 92(2):409-19

GRENGA PL, REALE G, COFONE C, MEDURI A, CERUTI P, GRENGA R
HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE
OPHTHAL PLAST RECONSTR SURG, 2009 MAR-APR; 25(2):123-5

CRESPI R, CAPPARÈ P, GHERLONE E
DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS
J PERIODONTOL, 2009 OCT; 80(10):1616-1621

RINNA C, REALE G, FORESTA E, MUSTAZZA MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4

CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N
PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88

SCARANO A, CARINCI F, ASSENZA B, PIATTELLI M, MURMURA G, PIATTELLI A
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES
CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30

PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58

SANTAGATA M, GUARINIELLO L, TARTARO G
A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N°:114-9

SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18

CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIUTTI E, GHERLONE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76

SLOTTE C, LINDFORS N, NANNMARK U
SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23

BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6

SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13

RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9

BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U
A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION
J PERIODONTOL, 2012 JUL;83(7):836-46

- ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33
- FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
 IMPLANT DENT, 2012 JUN;21(3):190-5
- BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
 CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
 INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
 EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (S), 103-107
- FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
 INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66
- FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RESSION TYPE DEFECTS: A CASE SERIES
 QUINTESSENCE INT, 2013;44(3):243-6
- SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
 INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
 INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- RODRIGUEZ JG, ELDBANY RM
VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION
 INT J ORAL MAXILLOFAC SURG, 2013 SEP;42(9):1060-6
- FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
 MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1;33(6):3506-13
- FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61
- TRAINI T, PIATTELLI A, CAPUTI S, DEGIDI M, MANGANO C, SCARANO A, PERROTTI V, IEZZI G
REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS
 CLIN IMPLANT DENT RELAT RES, 2015 FEB;17(1):150-62
- KOLMAS J, SZWAJA M, KOLODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENEIC BIOMATERIALS USED AS BONE SUBSTITUTES
 J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41
- PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARAUSSE C, BUTI J, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2013 WINTER;6(4):343-357
- FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M
PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS - AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):799-805
- BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- BARONE A, RICCI M, ROMANOS GE, TONELLI P, ALFONSI F, COVANI U
BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):823-30
- BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
 CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJIARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
 CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- GHEÑO E, PALERMO A, BUFFOLI B, RODELLA LF
THE EFFECTIVENESS OF THE USE OF XENOGENEIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES
 J OSSEOINTEGRATION 2014;6(2):37-42
- FALISI G, GALLI M, VITTORINI-VELASQUEZ P, GALLEGOS-RIVERA JC, MINASI R, DE BIASE A, DI PAOLO C
USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE
 ACTA ODONTOL LATINOAM 2013;26(3):167-172
- FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
 QUINTESSENCE INT, 2014 NOV-DEC; 45(10):853-60
- FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
 CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42
- SCARANO A, MURMURA G, SINJIARI B, ASSENZA B, SOLLAZZO V, SPINELLI G, CARINCI F
EXPANSION OF THE ALVEOLAR BONE CREST WITH ULTRASONIC SURGERY DEVICE: CLINICAL STUDY IN MANDIBLE
 INT J IMMUNOPATHOL PHARMACOL, 2011 APR-JUN; 24(2 SUPPL):71-5
- SCARANO A, PIATTELLI A, MURMURA G, IEZZI G, ASSENZA B, MANCINO C
DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGIC CASE SERIES
 INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9
- LORENZO G, BUTTARELLO GM, CHESSA G
CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION
 DENTISTRY, 2015, 5:2
- THALMAIR T, FICKL S, SCHNEIDER D, HINZE M, WACHTEL H
DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES - A VOLUMETRIC STUDY
 J CLIN PERIODONTOL, 2013 JUL;40(7):721-7
- MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATTELLI A, TRUBIANI O
OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON
 J PERIODONTAL RES, 2016 FEB;51(1):112-24
- SCARANO A, PIATTELLI A, IEZZI G, VARVARA G
SPONTANEOUS BONE FORMATION ON THE MAXILLARY SINUS FLOOR IN ASSOCIATION WITH SURGERY TO REMOVE A MIGRATED DENTAL IMPLANT: A CASE REPORT
 MINERVA STOMATOL, 2014 OCT;63(10):351-9
- BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL
 CLIN ORAL IMPLANTS RES, 2016 NOV;27(11):E105-E115
- ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- OZEL B, FINDIKIOGLU K, SEZGIN B, GUNAY K, BARUT I, OZMEN S
A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE
 J CRANIOMAXILLOFAC SURG, 2015 OCT;43(8):1583-8
- CORBELLA S, TASCHIERI S, WEINSTEIN R, DEL FABBRO M
HISTOMORPHOMETRIC OUTCOMES AFTER LATERAL SINUS FLOOR ELEVATION PROCEDURE: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS
 CLIN ORAL IMPLANTS RES, 2016 SEP;27(9):1106-22
- BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY
 CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112.
- CAKIR M, KARACA IR, AYŞEĞÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION
 CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85
- LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F
REGENERATION OF ATROPHIC CRESTAL RIDGES WITH RESORBABLE LAMINA: TECHNICAL NOTE
 J BIOL REGUL HOMEOST AGENTS 2015 JUL-SEP;29(3 SUPPL 1):97-100
- ETTORRE V, DE MARCO P, ZARA S, PERROTTI V, SCARANO A, DI CRESCENZO A, PETRINI M, HADAD C, BOSCO D, ZAVAN B, VALBONETTI L, SPOTO G, IEZZI G, PIATTELLI A, CATALDI A, FONTANA A
IN VITRO AND IN VIVO CHARACTERIZATION OF GRAPHENE OXIDE COATED PORCINE BONE GRANULES
 CARBON, JULY 2016, VOLUME 103, PAGES 291-298
- ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
 MINERVA STOMATOL, 2016 OCT;65(5):257-68
- SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
 BIOMED RES INT, 2016;2016:4086870
- FELICE P, ZUCHELLI G, CANNIZZARO G, BARAUSSE C, DIAZZI M, TRULLENQUE-ERIKSSON A, ESPOSITO M
IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 4-MONTH POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2016;9(3):233-247

ROMBOUTS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

BARONE A, TOTI P, MARCONCINI S, DERCHI G, MARCHIONNI S, COVANI U
ESTHETIC OUTCOME OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS BY CLINICIANS WITH OR WITHOUT EXPERIENCE: A MEDIUM-TERM RETROSPECTIVE EVALUATION
INT J ORAL MAXILLOFAC IMPLANTS, 2016;31(6)

TALLARICO M, KHANARI E, PISANO M, DE RIU G, TULLIO A, LORONI SM
SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):263-275

LOPEZ MA, MANZULLI N, CASALE M, ORMIANER Z, CARINCI F
THE USE OF RESORBABLE HETEROLOGOUS CORTICAL LAMINA AS A NEW SINUS LIFT FLOOR: A TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):75-79

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

ESPOSITO M, ZUCHELLI G, BARAUSSE C, PISTILLI R, TRULLENQUE-ERIKSSON A, FELICE P
FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 4-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(4):393-409

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U
CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTRE RANDOMIZED CONTROLLED CLINICAL TRIAL
J CLIN PERIODONTOL, 2017 FEB;44(2):204-214

ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P, BARONE A
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS
CURR PHARM BIOTECHNOL, 2017;18(1):64-75

BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

BARONE A, TOTI P, FUNEL N, CAMPANI D, COVANI U
EXPRESSION OF SP7, RUNX1, DLX5, AND CTNNB1 IN HUMAN MESENCHYMAL STEM CELLS CULTURED ON XENOGENIC BONE SUBSTITUTE AS COMPARED WITH MACHINED TITANIUM
IMPLANT DENT, 2014 AUG;23(4):407-15

ESPOSITO M, ZUCHELLI G, CANNIZZARO G, CHECCHI L, BARAUSSE C, TRULLENQUE-ERIKSSON, FELICE P
IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 1-YEAR POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(1):11-26

SCARANO A, CRINCOLI V, DI BENEDETTO A, COZZOLINO V, LORUSSO F, PODALIRI VULPIANI M, GRANO M, KALEMAJ Z, MORI G, GRASSI FR
BONE REGENERATION INDUCED BY BONE PORCINE BLOCK WITH BONE MARROW STROMAL STEM CELLS IN A MINIPIG MODEL OF MANDIBULAR "CRITICAL SIZE" DEFECT
STEM CELLS INT, 2017;2017:9082869

SCARANO A
TRADITIONAL POSTEXTRACTIVE IMPLANT SITE PREPARATION COMPARED WITH PRE-EXTRACTIVE INTERRADICULAR IMPLANT BED PREPARATION IN THE MANDIBULAR MOLAR REGION, USING AN ULTRASONIC DEVICE: A RANDOMIZED PILOT STUDY
INT J ORAL MAXILLOFAC IMPLANTS, 2017 MAY/JUN;32(3):655-660

BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLINICAL ORAL INVESTIGATIONS, 2018 JAN;22(1):505-513

IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44

FELICE P, BARAUSSE C, BARONE A, ZUCHELLI G, PIATTELLI M, PISTILLI R, IPPOLITO DR, SIMION M
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

IIDA T, CARNEIRO MARTINS NETO E, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP
INFLUENCE OF A COLLAGEN MEMBRANE POSITIONED SUBJACENT THE SINUS MUCOSA FOLLOWING THE ELEVATION OF THE MAXILLARY SINUS. A HISTOMORPHOMETRIC STUDY IN RABBITS
CLIN IMPLANT DENT RELAT RES, 2017 JUN 7

DE MARCO P, ZARA S, DE COLLI M, RADUNOVIC M, LAZOVIC V, ETORRE V, DI CRESCENZO A, PIATTELLI A, CATALDI A, FONTANA A
GRAPHENE OXIDE IMPROVES THE BIOCOMPATIBILITY OF COLLAGEN MEMBRANES IN AN IN VITRO MODEL OF HUMAN PRIMARY GINGIVAL FIBROBLASTS
BIOMED MATER, 2017 SEP 13;12(5):055005

MJIURITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR MATERIALS, 2017 SEP 8;10(9)

ROSSI R, FOCE E, SCOLAVINO S
THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT
J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41

CHECCHI V, FELICE P, ZUCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278

CRESPI R, CAPPARÈ P, GHERLONE E
COMPARISON OF MAGNESIUM-ENRICHED HYDROXYAPATITE AND PORCINE BONE IN HUMAN EXTRACTION SOCKET HEALING: A HISTOLOGIC AND HISTOMORPHOMETRIC EVALUATION
INT J ORAL MAXILLOFAC IMPLANTS, 2011 SEP-OCT;26(5):1057-62

CORBELLA S, TASCHIERI S, FRANCIETTI L, WEINSTEIN R, DEL FABBRO M
HISTOMORPHOMETRIC RESULTS AFTER POSTEXTRACTION SOCKET HEALING WITH DIFFERENT BIOMATERIALS: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2017 SEP/OCT;32(5):1001-1017

RADUNOVIC M, DE COLLI M, DE MARCO P, DI NISIO C, FONTANA A, PIATTELLI A, CATALDI A, ZARA S
GRAPHENE OXIDE ENRICHMENT OF COLLAGEN MEMBRANES IMPROVES DPSCS DIFFERENTIATION AND CONTROLS INFLAMMATION OCCURRENCE
J BIOMED MATER RES A, 2017 AUG;105(8):2312-2320

KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
BMC ORAL HEALTH, 2017 OCT 5;17(1):126

TROIANO G, ZHURAKIVSKA K, LO MUZIO L, LAINO L, CICCIO M, LO RUSSO L
COMBINATION OF BONE GRAFT AND RESORBABLE MEMBRANE FOR ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW, META-ANALYSIS AND TRIAL SEQUENTIAL ANALYSIS
J PERIODONTOL, 2017 SEP 12:1-17. EPUB AHEAD OF PRINT

ROSSI R, LONGO E, MIJIRITSKY E
A NEW INTERPRETATION OF GUIDED IMPLANT SURGERY TO ACHIEVE AN OPTIMAL RESULT IN THE ESTHETIC ZONES
MEDICAL RESEARCH ARCHIVES, 2017 APRIL, VOL. 5, ISSUE 4

BOLLE C, FELICE P, BARAUSSE C, PISTILLI V, TRULLENQUE-ERIKSSON A, ESPOSITO M
FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):31-47

ESPOSITO M, DAVÒ R, MARTI PAGES C, FERRER FUENTES A, BARAUSSE C, PISTILLI R, IPPOLITO DR, FELICE P
IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 4 MONTHS POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):11-28

GASTALDI G, FELICE P, PISTILLI V, BARAUSSE C, IPPOLITO DR, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 x 5 MM IMPLANTS WITH A NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):49-61

DIOMEDE F, D'AUROA M, GUGLIANDOLO A, MERCIARO I, ORSINI T, GATTA V, PIATTELLI A, TRUBIANI O, MAZZON E
BIOFUNCTIONALIZED SCAFFOLD IN BONE TISSUE REPAIR
INT J OF MOLECULAR SCIENCES, 2018, 19, 1022

SCARANO A, LORUSSO F, SANTOS DE OLIVEIRA P, MURMURA G, CARINCI F
DENTAL IMPLANTS DISPLACED INTO THE MANDIBULAR CORPUS: CLINICAL NOTE
JOURNAL OF DENTISTRY AND ORAL CARE, 2018, 4(1):8-12

FORABOSCO A, GHENO E, SPINATO S, GARUTI G, FORABOSCO E, CONSOLO U
CONCENTRATED GROWTH FACTORS IN MAXILLARY SINUS FLOOR AUGMENTATION: A PRELIMINARY CLINICAL COMPARATIVE EVALUATION
INT J OF GROWTH FACTORS AND STEM CELLS IN DENTISTRY, 2018;1:2-7

PISTILLI R, FELICE P, CANNIZZARO G, PIATTELLI M, CORVINO V, BARAUSSE C, BUTI J, SOARDI E, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG, 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR POST-LOADING RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2013;6(4):359-372

- FELICE P, BARAUSSE C, PISTILLI V, PIATTELLI M, IPPOLITO DR, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR POST-LOADING RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(2):175-187
- DAVÓ R, FELICE P, PISTILLI R, BARAUSSE C, MARTI PAGES C, FERRER FUERTES A, IPPOLITO DR, ESPOSITO M
IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2018;11(1):145-161
- OMORI Y, SILVA ER, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP
REPOSITION OF THE BONE PLATE OVER THE ANTROSTOMY IN MAXILLARY SINUS AUGMENTATION: A HISTOMORPHOMETRIC STUDY IN RABBITS
CLIN ORAL IMPLANTS RES, 2018 AUG;29(8):821-834
- CANULLO L, GENOVA T, NAENNIC N, NAKAJIMA Y, MASUDA K, MUSSANO F
PLASMA OF ARGON ENHANCES THE ADHESION OF MURINE OSTEOBLASTS ON DIFFERENT GRAFT MATERIALS
ANN ANAT, 2018 JUL;218:265-270
- RAMÍREZ FERNÁNDEZ MP, MAZÓN P, GEHRKE SA, CALVO GUIRADO JL, DE AZA PN
COMPARISON OF TWO XENOGRFT MATERIALS USED IN SINUS LIFT PROCEDURES: MATERIAL CHARACTERIZATION AND IN VIVO BEHAVIOR
MATERIALS, 2017 JUN 7;10(6)
- DIOMEDE F, ZINI N, GATTA V, FULLE S, MERCIARO I, D'AUROA M, LA ROVERE RML, TRAINI T, PIZZICANNELLA J, BALLERINI P, PIATTELLI A, TRUBIANI O
HUMAN PERIODONTAL LIGAMENT STEM CELLS CULTURED ONTO CORTICO-CANCELLOUS SCAFFOLD DRIVE BONE REGENERATIVE PROCESS
EUR CELL MATER, 2016 SEP 16;32:181-201
- DEVELIOGLU H, OZCAN G, GULTEKIN SE, SENGUVEN B, YILDIRIM A
THE SHORT-TERM EFFECTS OF VARIOUS XENOGRFTS ON BONE HEALING IN RATS CRANIAL DEFECTS
BIOMEDICAL RESEARCH, 2018;29(8): 1598-1602
- BRUNELLI G, CARINCI F, GIRARDI A, PALMIERI A, CACCIANIGA GL, SOLLAZZO V
OSTEOBIOL® EFFECT ON DENTAL PULP DERIVED STEM CELLS EUR J INFLAMMATION, VOL. 10, NO. 1 (S), 27-30
- MAZZONI S, MOHAMMADI S, TROMBA G, DIOMEDE F, PIATTELLI A, TRUBIANI O, GIULIANI A
ROLE OF CORTICO-CANCELLOUS HETEROLOGOUS BONE IN HUMAN PERIODONTAL LIGAMENT STEM CELLS XENO-FREE CULTURE STUDIED BY SYNCHROTRON RADIATION PHASE-CONTRAST MICROTOMOGRAPHY
INT J MOL SCI, 2017;18;364
- LAURITANO D, CARINCI F, ZOLLINO I, HASSANIPOUR A, SAGGESE V, PALMIERI A, GIRARDI A, CURA F, PIRAS A, ZAMBONI P, BRUNELLI G
OSTEOBIOL® ENHANCES OSTEOGENIC DIFFERENTIATION IN BONE MARROW DERIVED STEM CELLS
EUR J INFLAMMATION, VOL. 10, NO. 1 (S3), 83-88
- MATÉ SANCHEZ DE VAL J, MAZON P, PIATTELLI A, CALVO GUIRADO J, MAREQUE BUENO J, GRANERO MARIN J, DE AZA P
COMPARISON AMONG THE PHYSICAL PROPERTIES OF CALCIUM PHOSPHATE-BASED BONE SUBSTITUTES OF NATURAL OR SYNTHETIC ORIGIN
INT J APPLIED CERAMIC TECHNOL, 2018;15(4):930-937
- NAKAJIMA Y, PIATTELLI A, IEZZI G, FORTICH MESA N, FERRI M, BOTTICELLI D
INFLUENCE OF THE PRESENCE OF ALVEOLAR MUCOSA AT IMPLANTS: A HISTOLOGICAL STUDY IN HUMANS
IMPLANT DENT, 2018;27(2):193-201
- CASCONI P, VELLONE V, RAMIERI V, BASILE E, TARSITANO A, MARCHETTI C
RECONSTRUCTION OF THE ADULT HEMIFACIAL MICROSNOMIA PATIENT WITH TEMPOROMANDIBULAR JOINT TOTAL JOINT PROSTHESIS AND ORTHOGNATIC SURGERY
CASE REP SURG, 2018;1-6
- AMR AEH, ABDEL GHAFFAR KA, ABUEL-ELA HA, ABD ELHAMID ES
XENOGENIC FLEXIBLE BONE LAMINA GRAFT: A SUCCESSFUL ALTERNATIVE TO THE AUTOGENOUS ONLY BONE BLOCK GRAFT IN ALVEOLAR RIDGE AUGMENTATION: A CLINICAL, RADIOGRAPHIC AND HISTOLOGICAL EVALUATION
JOURNAL OF DENTAL TREATMENT AND ORAL CARE, 2017;1(1);104
- ATTIA AM
CLINICAL AND RADIOGRAPHIC EVALUATION OF CORTICO-CANCELLOUS BONE MIX XENOGRFT (OSTEOBIOL GEN-OS) IN THE TREATMENT OF HUMAN PERIODONTAL INTRABONY DEFECTS
EGYPTIAN DENTAL JOURNAL, 2017;63;1-10
- BECHARA S, KUBILIUS R, VERONESI G, PIRES JT, SHIBLI JA, MANGANO FG
SHORT (6-MM) DENTAL IMPLANTS VERSUS SINUS FLOOR ELEVATION AND PLACEMENT OF LONGER (>10-MM) DENTAL IMPLANTS: A RANDOMIZED CONTROLLED TRIAL WITH A 3-YEAR FOLLOW-UP
CLIN ORAL IMPLANTS RES, 2017;28;1097-1107
- CHANDRASEKARAN B, SURESH N, MUTHUSAMY S
PLATELET-RICH FIBRIN WITH BONE GRAFTS FOR REGENERATION OF BONY DEFECT FOLLOWING EXTRACTION OF SUPERNUMERARY TEETH: A CASE REPORT
CHIN J DENT RES, 2017;20(4):231-234
- CHIRILA L, ROTARU C, FILIPOV I, SANDULESCU M
MANAGEMENT OF ACUTE MAXILLARY SINUSITIS AFTER SINUS BONE GRAFTING PROCEDURES WITH SIMULTANEOUS DENTAL IMPLANTS PLACEMENT - A RETROSPECTIVE STUDY
BMC INFECT DIS, 2016 MAR 8;16 SUPPL 1:94
- EKSTEIN J, TANDELICH M, NART J, CALVO GUIRADO JL, SHAPIRA L
MARGINAL BONE LEVEL AROUND CONICAL CONNECTION TAPERED IMPLANTS WITH PLATFORM SWITCHING: A MULTICENTER RETROSPECTIVE STUDY AT 14 MONTHS FOLLOW-UP
J OSSEOINTEGRATION, 2016 JAN-APR;8(1)
- DEL CORSO M, DOHAN EHRENFEST DM
IMMEDIATE IMPLANTATION AND PERI-IMPLANT NATURAL BONE REGENERATION (NBR) IN THE SEVERELY RESORBED POSTERIOR MANDIBLE USING LEUKOCYTE- AND PLATELET-RICH FIBRIN (L-PRF): A 4-YEAR FOLLOW-UP
POSEIDO, 2013;1(2);85-92
- NOAMI SA, ELMOSY K, ASKAR N
EVALUATION OF PRE-HYDRATED COLLAGENATED CORTICO-CANCELLOUS GRANULES (MP3®) IN AUGMENTATION OF THE MAXILLARY SINUS (PRELIMINARY STUDY)
J OF DENT, ORAL DISORD AND THERAPY, 2014;2(3): 1-8
- BARONE A, TOTI P, PIATTELLI A, IEZZI G, DERCHI G, COVANI U
EXTRACTION SOCKET HEALING IN HUMANS AFTER RIDGE PRESERVATION TECHNIQUES: COMPARISON BETWEEN FLAPLESS AND FLAPPED PROCEDURES IN A RANDOMIZED CLINICAL TRIAL
J PERIODONTOL, 2014 JAN;85(1):14-23
- LOPEZ MA, LICO S, CASALE M, ORMANIER Z, CARINCI F
THE USE OF VARIOUS BIOMATERIALS IN COMPUTER-GUIDED CRESTAL SINUS LIFT PROCEDURES. A REPORT ON TWO CASE STUDIES WITH VOLUME COMPARISON
ORAL IMPLANTOL, 2016;APR-JUN 9(2)
- MEHL C, BOSCH T
IMPLANTOLOGY AND PROTHODONTICS AT CROSSROADS - CLASSIC VS MODERN TREATMENT CONCEPTS
JSM DENTAL SURGERY, 2016;1(1);1003
- NEMTOI A, DANILA V, DRAGAN E, PASCA S, NEMTOI A, CONSTANTIN M, SAVA A, HABA D
THE EFFECTS OF INSULIN AND STRONTIUM RANELATE ON GUIDED BONE REGENERATION IN DIABETIC RATS
REVISTA DE CHIMIE, 2017;68(4);693-697
- RAMIREZ FERNANDEZ MP, GEHRKE SA, ALBACETE MARTINEZ CP, CALVO GUIRADO JL, DE AZA PN
SEM-EDX STUDY OF THE DEGRADATION PROCESS OF TWO XENOGRFT MATERIALS USED IN SINUS LIFT PROCEDURES
MATERIALS, 2017;10;542
- ZITA GOMES R, DE VASCONCELOS MR, LOPES GUERRA IM, DE ALMEIDA RAB, DE CAMPOS FELINO AC
IMPLANT STABILITY IN THE POSTERIOR MAXILLA: A CONTROLLED CLINICAL TRIAL
BIOMED RES INT, 2017;6825213
- COVANI U, CANULLO L, TOTI P, ALFONSI F, BARONE A
TISSUE STABILITY OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: A 5-YEAR PROSPECTIVE SINGLE-COHORT STUDY
J PERIODONTOL, 2014 SEP;85(9):E323-32
- DIOMEDE F, D'AUROA M, GUGLIANDOLO A, MERCIARO I, ETORRE V, BRAMANTI A, PIATTELLI A, GATTA V, MAZZON E, FONTANA A, TRUBIANI O
A NOVEL ROLE IN SKELETAL SEGMENT REGENERATION OF EXTRACELLULAR VESICLES RELEASED FROM PERIODONTAL-LIGAMENT STEM CELLS
INT J NANOMEDICINE, 2018 JUN 29;13:3805-3825
- KAWAKAMI S, LANG NP, IIDA T, FERRI M, APAZA ALCCAYHUAMAN KA, BOTTICELLI D
INFLUENCE OF THE POSITION OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED WITH CONE-BEAM COMPUTED TOMOGRAPHY: A RANDOMIZED CLINICAL TRIAL
JOURNAL INVESTIG CLIN DENT, 2018 NOV;9(4):E12362
- SCARANO A, DE OLIVEIRA PS, TRAINI T, LORUSSO F
SINUS MEMBRANE ELEVATION WITH HETEROLOGOUS CORTICAL LAMINA: A RANDOMIZED STUDY OF A NEW SURGICAL TECHNIQUE FOR MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT
MATERIALS, 2018 AUG 17;11(8)
- FIGLIUZZI MM, GIUDICE A, CRISTOFARO MG, PACIFICO D, BIAMONTE P, FORTUNATO L
POSTEXTRACTIVE IMPLANTS IN AESTHETIC AREAS: EVALUATION OF PERIMPLANT BONE REMODELING OVER TIME
ANN STOMATOL, 2015 MAY 18;6(1):29-34
- KIVOVICS M, SZABÓ BT, NÉMETH O, TARI N, DORI F, NAGY P, DOBÓ-NAGY C, SZABÓ G
MICROARCHITECTURAL STUDY OF THE AUGMENTED BONE FOLLOWING RIDGE PRESERVATION WITH A PORCINE XENOGRFT AND A COLLAGEN MEMBRANE: PRELIMINARY REPORT OF A PROSPECTIVE CLINICAL, HISTOLOGICAL, AND MICRO-COMPUTED TOMOGRAPHY ANALYSIS
INT J ORAL MAXILLOFAC SURG, 2017 FEB;46(2):250-260
- MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A
CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL
CLIN IMPL DENT RELAT RES, 2018 DEC;20(6):906-914
- MARCONCINI S, COVANI U, GIAMMARINARO E, VELASCO-ORTEGA E, DE SANTIS D, ALFONSI F, BARONE A
CLINICAL SUCCESS OF DENTAL IMPLANTS PLACED IN POSTERIOR MANDIBLE AUGMENTED WITH INTERPOSITIONAL BLOCK GRAFT: 3-YEAR RESULTS FROM A PROSPECTIVE COHORT CLINICAL STUDY
J ORAL MAXILLOFAC SURG. 2019 FEB;77(2):289-298
- IIDA T, RICARDO SILVA E, LANG NP, APAZA ALCCAYHUAMAN KA, DANIELE BOTTICELLI D, XAVIER SP
HISTOLOGICAL AND MICRO-COMPUTED TOMOGRAPHY EVALUATIONS OF NEWLY FORMED BONE AFTER MAXILLARY SINUS AUGMENTATION USING A XENOGRFT WITH SIMILAR DENSITY AND MINERAL CONTENT OF BONE: AN EXPERIMENTAL STUDY IN RABBITS
CLIN EXP DENT RES, 2018;1-7
- BERNARDI S, GATTO R, SEVERINO M, BOTTICELLI G, CARUSO S, RASTELLI C, LUPI E, ROIAS AQ, IACOMINO E, FALISI G
SHORT VERSUS LONGER IMPLANTS IN MANDIBULAR ALVEOLAR RIDGE AUGMENTED USING OSTEOGENIC DISTRACTION: ONE-YEAR FOLLOW-UP OF A RANDOMIZED SPLIT-MOUTH TRIAL
J ORAL IMPLANTOL, 2018 JUN;44(3):184-191
- DIKER N, SARICA H, CUMBUL A, KILIC E
EFFECTS OF SYSTEMIC ERYTHROPOIETIN TREATMENT AND HETEROGENEOUS XENOGRFT IN COMBINATION ON BONE REGENERATION OF A CRITICAL-SIZE DEFECT IN AN EXPERIMENTAL MODEL
J CRANIOMAXILLOFAC SURG, 2018 NOV;46(11):1919-1923
- SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F
A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592
- SENESE O, BOUTREMAN E, GOSSIAUX C, LOEB I, DEQUANTER D
RETROSPECTIVE ANALYSIS OF 79 PATIENTS WITH ORBITAL FLOOR FRACTURE: OUTCOMES AND PATIENT-REPORTED SATISFACTION
ARCH CRANIOFAC SURG, 2018 JUN;19(2):108-113

GENOVA T, PESCE P, MUSSANO F, TANAKA K, CANULLO L
**THE INFLUENCE OF BONE-GRAFT
BIO-FUNCTIONALIZATION WITH PLASMA OF ARGON
ON BACTERIAL CONTAMINATION**
J BIOMED MATER RES A, 2019 JAN;107(1):67-70

ESPOSITO M, BARAUSSE C, PISTILLI R, PIATTELLI M, DI SIMONE S,
IPPOLITO DR, FELICE P
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH
PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A
NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED
TITANIUM SURFACE OR BY LONGER IMPLANTS IN
AUGMENTED BONE. FIVE-YEAR RESULTS FROM A
RANDOMISED CONTROLLED TRIAL**
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

FELICE P, PISTILLI R, BARAUSSE C, PIATTELLI M, BUTI J,
ESPOSITO M
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH
PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE
IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED
BONE. FIVE-YEAR POST-LOADING RESULTS FROM A
WITHIN-PERSON RANDOMISED CONTROLLED TRIAL**
INT J OF ORAL IMPLANTOL, 2019;12(1):57-72

DI CARLO R, ZARA S, VENTRELLA A, SIANI G, DA ROS T, IEZZI G,
CATALDI A, FONTANA A
**COVALENT DECORATION OF CORTICAL MEMBRANES
WITH GRAPHENE OXIDE AS A SUBSTRATE FOR DENTAL
PULP STEM CELLS**
NANOMATERIALS, 2019;9:604

CABALLÉ-SERRANO J, MUNAR-FRAUJA A, DELGADO L, PÉREZ R,
HERNÁNDEZ-ALFARO F
**PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER
MEMBRANES FOR BONE REGENERATION**
J MECH BEHAV BIOMED, 2019;97:13-20

MATOH U, PETELIN M, GASPERIS R
**SPLIT-MOUTH COMPARISON OF CORONALLY
ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR
COLLAGEN MATRIX FOR TREATMENT OF ISOLATED
GINGIVAL RECESSIONS**
INT J PERIODONTICS RESTORATIVE DENT,
2019;39(3):439-446

AMBROZEWICZ E, MUSZYNSKA M, TOKAJUK G,
GRYNKIEWICZ G, ŻARKOWICZ N, SKRZYDLEWSKA
**BENEFICIAL EFFECTS OF VITAMINS K AND D3 ON REDOX
BALANCE OF HUMAN OSTEOBLASTS CULTURED WITH
HYDROXYAPATITE-BASED BIOMATERIALS
CELLS, 2019;8:325**

POLIS-YANES C, CADENAS-SEBASTIÁN C, GUALVAQUÉS P,
AYUSO-MONTERO R, MARÍ-ROIG A, LÓPEZ-LÓPEZ J
**GUIDED BONE REGENERATION OF AN ATROPHIC
MAXILLA USING HETEROLOGOUS CORTICAL LAMINA
CASE REP DENT, 2019; 5216362**

KAWAKAMI S, LANG NP, FERRI M, APAZA ALCCAYHUAMAN KA,
BOTTICELLI D
**INFLUENCE OF THE HEIGHT OF THE ANTROSTOMY IN SINUS
FLOOR ELEVATION ASSESSED BY CONE BEAM COMPUTED
TOMOGRAPHY- A RANDOMIZED CLINICAL TRIAL**
INT J ORAL MAXILLOFAC IMPLANTS, 2019;34(1):223-232

CHECCHI V, GASPARRO R, PISTILLI R, CANULLO L, FELICE P
**CLINICAL CLASSIFICATION OF BONE AUGMENTATION
PROCEDURE FAILURES IN THE ATROPHIC ANTERIOR
MAXILLAE: ESTHETIC CONSEQUENCES AND TREATMENT
OPTIONS**
BIOMED RES INT, 2019;4386709

FISCHER K R, TESTORI T, WACHTEL H, MÜHLEMANN S, HAPPE A,
DEL FABBRO M
**SOFT TISSUE AUGMENTATION APPLYING A
COLLAGENATED PORCINE DERMAL MATRIX DURING
SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER
CASE SERIES**
CLIN IMPLANT DENT RELAT RES.,2019;1-8

HIROTA A, LANG NP, FERRI M, FORTICH MESA N,
APAZA ALCCAYHUAMAN KA, BOTTICELLI D
**TOMOGRAPHIC EVALUATION OF THE INFLUENCE OF
THE PLACEMENT OF A COLLAGEN MEMBRANE
SUBJACENT TO THE SINUS MUCOSA DURING
MAXILLARY SINUS FLOOR AUGMENTATION: A
RANDOMIZED CLINICAL TRIAL**
INT J IMPLANT DENT.,2019; 5(1):31

TANAKA K, IEZZI G, PIATTELLI A, FERRI M, MESA NF,
APAZA ALCCAYHUAMAN KA, BOTTICELLI D
**SINUS FLOOR ELEVATION AND ANTROSTOMY HEALING:
A HISTOMORPHOMETRIC CLINICAL STUDY IN HUMANS**
IMPLANT DENT., 2019; AUG 14

JEANNEAU C, LE FOURNIS C, ABOUT I
**XENOGENIC BONE FILLING MATERIALS MODULATE
MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE
COMPLEMENT C5A**
CLIN. ORAL INVESTIG.; 2019 OCT 23

RAMANAUSKAITE A, BORGES T, ALMEIDA BL, CORREIA A
**DENTAL IMPLANT OUTCOMES IN GRAFTED SOCKETS: A
SYSTEMATIC REVIEW AND META-ANALYSIS**
J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E8

FARIA-ALMEIDA R, ASTRAMSKAITE-JANUSEVICIENE I, PUISYS A,
CORREIA F
**EXTRACTION SOCKET PRESERVATION WITH OR
WITHOUT MEMBRANES, SOFT TISSUE INFLUENCE ON
POST EXTRACTION ALVEOLAR RIDGE PRESERVATION: A
SYSTEMATIC REVIEW**
J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E5

ADILGÖLU S, GIRAY CB, KULAC I, USUBUTUN A, AKTAS A
**CLINICAL AND HISTOPATHOLOGICAL COMPARATIVE
STUDY OF TWO EQUINE-DERIVED BONE GRAFT: A
HUMAN STUDY**
J PAK MED ASSOC., 2019; 69: 1617

ROSSI R, FOCE E
**RECONSTRUCTION OF A HORIZONTAL AND VERTICAL
BONE DEFECT USING THE CORTICAL LAMINA
TECHNIQUE**
MED. RES. ARCH, 2019; VOL 7 (11)

CANULLO L, GENOVA T, RAKIC M, SCULEAN A, MIRON R,
MUZZI M, CAROSSA S, MUSSANO F
**EFFECTS OF ARGON PLASMA TREATMENT ON THE
OSTEOCONDUCTIVITY OF BONE GRAFTING MATERIALS**
CLIN ORAL INVESTIG, 2020 AUG;24(8):2611-2623

ROSSI R, CONTI A, BERTAZZO D, PILLONI A
**HORIZONTAL RIDGE AUGMENTATION WITH THE
CORTICAL LAMINA TECHNIQUE: A CASE REPORT**
MOD RES DENT; 2019; 4(4)

FELICE P, BARAUSSE C, BUTI J, GESSAROLI M, ESPOSITO M
**IMMEDIATE, EARLY (6 WEEKS) AND DELAYED (4
MONTHS) SINGLE POST-EXTRACTIVE IMPLANTS: 3-YEAR
POST-LOADING DATA FROM A RANDOMISED
CONTROLLED TRIAL**
CLINICAL TRIALS IN DENTISTRY; 2019;01(1):5-23

BARAUSSE C, FELICE P, PISTILLI R, BUTI J, ESPOSITO M
**POSTERIOR JAW REHABILITATION USING PARTIAL
PROSTHESES SUPPORTED BY IMPLANTS 4.0 X 4.0 MM OR
LONGER: THREE-YEAR POST-LOADING RESULTS OF A
MULTICENTRE RANDOMISED CONTROLLED TRIAL**
CLINICAL TRIALS IN DENTISTRY; 2019;01(1):25-36

VERARDI S, ORSINI M, LOMBARDI T, AUSENDA F, TESTORI T,
PULICI A, OREGLIA F, VALENTE NA, STACCHI C
**COMPARISON BETWEEN TWO DIFFERENT TECHNIQUES
FOR PERI-IMPLANT SOFT TISSUE AUGMENTATION:
PORCINE DERMAL MATRIX GRAFT VERSUS TENTING
SCREW**
J PERIODONTOL. 2020 AUG;91(8):1011-1017

KIZILASLAN S, KARABUDA ZC, ÖLGAC V
**THE EFFECT OF CONCENTRATED GROWTH FACTOR ON
CALVARIAL BONE IN DIABETIC HEALING**
J CRANIOFAC SURG 2020 JAN/FEB;31(1):158-161

ASLAN S, BUDULENI N, CORTELLINI P
**ENTIRE PAPILLA PRESERVATION TECHNIQUE IN THE
REGENERATIVE TREATMENT OF DEEP INTRABONY
DEFECTS: 1-YEAR RESULTS**
J CLIN PERIODONTOL; 2017 SEP; 44(9):926-932

FAVERO G, VIÑA-ALMUNIA J, CARDA C, MARTÍN DE LLANO JJ,
GARCÍA-MIRA B, SOTO-PEÑALOZA D, PEÑARROCHA-DIAGO M,
BOTTICELLI D
**INFLUENCE OF THE USE OF AUTOGENOUS BONE
PARTICLES TO CLOSE THE ACCESS WINDOW AFTER
MAXILLARY SINUS FLOOR AUGMENTATION: AN
EXPERIMENTAL STUDY IN RABBITS**
INT J IMPLANT DENT. 2020 MAR 4;6(1):9

TALLARICO M, KHANARI E, PISANO M, GATTI F, MELONI SM
**MOLAR REPLACEMENT WITH 7 MM-WIDE DIAMETER
IMPLANTS: TO PLACE THE IMPLANT IMMEDIATELY OR TO
WAIT 4 MONTHS AFTER SOCKET PRESERVATION? 1 YEAR
AFTER LOADING RESULTS FROM A RANDOMISED
CONTROLLED TRIAL**
EUR J ORAL IMPLANTOL, 2017;10(2):169-178

FELICE P, BARAUSSE C, PISTILLI R, KALEMAJ Z, ESPOSITO M
**FOUR-MM-LONG VERSUS LONGER IMPLANTS IN
AUGMENTED BONE IN ATROPHIC POSTERIOR JAWS:
THREE-YEAR POST-LOADING RESULTS FROM A
MULTICENTRE RANDOMISED CONTROLLED TRIAL**
CLINICAL TRIALS IN DENTISTRY, 2020;02(1):5-26

GIULIANI A, SENA G, TROMBA G, MAZZON E, FONTANA A,
DIOMEDE F, PIATTELLI A, TRUBIANI O
**COULD THE ENRICHMENT OF A BIOMATERIAL WITH
CONDITIONED MEDIUM OR EXTRACELLULAR VESICLES
MODIFY BONE-REMODELING KINETICS DURING A
DEFECT HEALING? EVALUATIONS ON RAT CALVARIA
WITH SYNCHROTRON-BASED MICROTOMOGRAPHY**
APPL SCI, 2020, 10, 2336

BALDI N, BUTI J, MENSÌ M, ALFONSI F, CINQUINI C, TONELLI P,
BARONE A
**XENOGENIC DERMAL MATRIX VERSUS AUTOLOGOUS
CONNECTIVE TISSUE GRAFT VERSUS NO GRAFT AT
ABUTMENT CONNECTION FOR IMPROVING
AESTHETICS: 6-MONTH OUTCOMES OF A RANDOMISED
CONTROLLED TRIAL**
CLINICAL TRIALS IN DENTISTRY, 2020;02(2):49-62

FISCHER KR, GÖTZ W, KAUFFMANN F, SCHMIDLIN PR,
FRIEDMANN A
**RIDGE PRESERVATION OF COMPROMISED EXTRACTION
SOCKETS APPLYING A SOFT CORTICAL MEMBRANE: A
CANINE PROOF-OF-PRINCIPLE EVALUATION**
ANN ANAT, 2020 SEP;231:151524

TOLEDANO M, ASADY S, TOLEDANO OSORIO M, GARCÍA
GODOY F, SERRERA FIGALLO MA, BENÍTEZ GARCÍA JA,
OSORIO R
**DIFFERENTIAL BIODEGRADATION KINETICS OF
COLLAGEN MEMBRANES FOR BONE REGENERATION**
POLYMERS, 2020 JUN 4;12(6):1290

LUONGO R, SGARAMELLA N, TRAINI T, BUGECA C
**GRAFTLESS MAXILLARY SINUS FLOOR AUGMENTATION
WITH SIMULTANEOUS PORCINE BONE LAYER
INSERTION: A 1- TO 5-YEAR FOLLOW-UP STUDY**
INT J ORAL MAXILLOFAC IMPLANTS, JUL/AUG
2020;35(4):808-815

SCARANO A, CICCARESE S, AMUSO D, MORTELLARO C,
LORUSSO F
**CORTICAL BONE LAMINA APPROACH FOR MANDIBULAR
LARGE CYSTIC DEFECT: A CASE REPORT**
J BIOL REGUL HOMEOST AGENTS, NOV-DEC 2019;33(6
SUPPL. 2):85-91

ASLAN S, BUDUNELI N, CORTELLINI P
**ENTIRE PAPILLA PRESERVATION TECHNIQUE: A NOVEL
SURGICAL APPROACH FOR REGENERATIVE TREATMENT
OF DEEP AND WIDE INTRABONY DEFECTS**
INT J PERIODONTICS RESTORATIVE DENT, MAR/APR
2017;37(2):227-233

FELICE P, BARAUSSE C, DAVÓ R, PISTILLI R, MARTI PAGES C,
FERRER FUERTES A, FERRI A
**IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VERSUS
CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED
ATROPHIC MAXILLAE: THREE-YEAR POST-LOADING
RESULTS FROM A MULTICENTRE RANDOMISED
CONTROLLED TRIAL**
CLINICAL TRIALS IN DENTISTRY, 2020;02(3):5-25

ESPOSITO M, BARAUSSE C, PISTILLI R, BELLINI P, BUTI J, FELICE P
**IMMEDIATE LOADING OF 3 MM-DIAMETER IMPLANTS AS
AN ALTERNATIVE HORIZONTAL BONE AUGMENTATION
FOR PLACING 4 MM-DIAMETER IMPLANTS: ONE-YEAR
POST-LOADING RESULTS FROM A MULTICENTRE
RANDOMISED CONTROLLED TRIAL**
CLINICAL TRIALS IN DENTISTRY, 2020;02(3):61-76

BERGMANN M, JEANNEAU C, GIRAUD T, RICHARD G, ABOUT I
**COMPLEMENT ACTIVATION LINKS INFLAMMATION TO
DENTAL TISSUE REGENERATION**
CLIN ORAL INVESTIG, 2020 DEC;24(12):4185-4196

ROSSI R, GHEZZI C, TOMECEK M
**CORTICAL LAMINA: A NEW DEVICE FOR THE TREATMENT
OF MODERATE AND SEVERE TRIDIMENSIONAL BONE
AND SOFT TISSUE DEFECTS**
INT J ESTHET DENT, 2020;15(4):454-473

IEZZI G, PERROTTI V, PIATTELLI A, D'ARCANGELO C, ALTIERI F,
CASSETTA M
**SPLIT CREST TECHNIQUE FOR IMPLANT TREATMENT OF
AGENESIS OF THE UPPER LATERAL INCISORS: RESULTS
OF A RANDOMIZED PILOT HISTOLOGICAL AND
CLINICAL STUDY AT 24-MONTH FOLLOW-UP**
BRAZ ORAL RES, 2020 SEP 4;34:E118

ESPOSITO M, BARAUSSE C, BONIFAZI L, PIATTELLI M, PISTILLI R, FERRI A, FELICE P
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VERSUS SOCKET PRESERVATION AND DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN THE MOLAR REGION: 5-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
CLINICAL TRIALS IN DENTISTRY, 2020;02(4):21-37

ARAGONESES J, SUÁREZ A, RODRÍGUEZ C, ARAGONESES JM
CLINICAL AND HISTOLOGICAL DIFFERENCES BETWEEN GUIDED TISSUE REGENERATION WITH ACELLULAR DERMAL MATRIX OF PORCINE ORIGIN AND AUTOLOGOUS CONNECTIVE TISSUE: AN ANIMAL STUDY
MATERIALS, 2021 JAN 7;14(2):272

TH ELASKARY A, GAWEESH YY, EL TANTAWI M, MAEBED MA
VESTIBULAR SOCKET THERAPY: A NOVEL APPROACH FOR IMPLANT PLACEMENT IN DEFECTIVE FRESH EXTRACTION SOCKETS WITH OR WITHOUT ACTIVE SOCKET INFECTION (ONE-ARM COHORT STUDY)
INT J ORAL MAXILLOFAC IMPLANTS, JAN-FEB 2021;36(1):146-153

FALACHO RI, PALMA PJ, MARQUES JA, FIGUEIREDO MH, CAMELO F, DIAS I, VIEGAS C, GUERRA F
COLLAGENATED PORCINE HETEROLOGOUS BONE GRAFTS: HISTOMORPHOMETRIC EVALUATION OF BONE FORMATION USING DIFFERENT PHYSICAL FORMS IN A RABBIT CANCELLOUS BONE MODEL
MOLECULES, 2021 MAR 2;26(5):1339

IGLESIAS VELÁZQUEZ Ó, TRESGUERRES FGF, LECO BERROCAL I, TRESGUERRES IF, LÓPEZ PINTOR RM, CARBADILLO J, LÓPEZ QUILES J, TORRES J
SPLIT BONE BLOCK TECHNIQUE: 4-MONTH RESULTS OF A RANDOMISED CLINICAL TRIAL COMPARING CLINICAL AND RADIOGRAPHIC OUTCOMES BETWEEN AUTOGENOUS AND XENOGENIC CORTICAL PLATES
INT J OF ORAL IMPLANTOLOGY, 2021;14(1):41-52

FOCE E, ROSSI R
HORIZONTAL BONE AUGMENTATION WITH THE CORTICAL LAMINA TECHNIQUE: A CASE REPORT
BIOMED J SCI & TECH RES, 35(3)-2021

LOPEZ MA, PASSARELLI PC, RELLA E, ALTAMURA FR, SANTACROCE L, CASALE M, D'ADDONA A
SINUS PACK FOR MAXILLARY SINUS AUGMENTATION: A NEW TECHNIQUE
J OSSEOINTEGR, 2021;13(2):1

ELASKARY A, MAEBED M, ABD-ELWAHAB RAD I
VESTIBULAR SOCKET THERAPY WITH IMMEDIATE IMPLANT PLACEMENT FOR MANAGING COMPROMISED FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE-ARM CLINICAL STUDY
INT J ORAL IMPLANTOL (BERL), 2021 AUG 20;14(3):307-320

SCHUH PL, WACHTEL H, BEUER F, GOKER F, DEL FABBRIO M, FRANCIOTTI L, TESTORI T
MULTI-LAYER TECHNIQUE (MLT) WITH PORCINE COLLAGENATED CORTICAL BONE LAMINA FOR BONE REGENERATION PROCEDURES AND IMMEDIATE POST-EXTRACTION IMPLANTATION IN THE ESTHETIC AREA: A RETROSPECTIVE CASE SERIES WITH A MEAN FOLLOW-UP OF 5 YEARS
MATERIALS, 2021 SEP 9;14(18):5180

ORTIZ PUIGPELAT O, SIMÕES A, CABALLÉ SERRANO J, HERNÁNDEZ ALFARO F
BLOOD ABSORPTION CAPACITY OF DIFFERENT XENOGRAFT BONE SUBSTITUTES: AN IN-VITRO STUDY
J CLIN EXP DENT, 2019 NOV 1;11(11):E1018-E1024

CRESPI R, TOTI P, CRESPI G, COVANI U, BREVI B, MENCHINI FABRIS GB
BONE REMODELING AROUND IMPLANTS PLACED AFTER SOCKET PRESERVATION: A 10-YEAR RETROSPECTIVE RADIOLOGICAL STUDY
INT J IMPLANT DENT, 2021 JUL 29;7(1):74

SCALA A, VIÑA ALMUNIA J, CARCA C, MARTÍN DE LLANO JJ, SOTO PEÑALOZA D, PEÑARROCHA DIAGO M, PEÑARROCHA DIAGO M, BOTTICELLI D
SEQUENTIAL HEALING OF THE ELEVATED SINUS FLOOR WITH DIFFERENT SIZE OF ANTROSTOMY: A HISTOMORPHOMETRIC STUDY IN RABBITS
ORAL MAXILLOFAC SURG, 2020 DEC;24(4):403-410

BERNARDELLO F, LOMBARDI T, STACCHI C
CLEARANCE OF BONE SUBSTITUTE IN GEL FORM ACCIDENTALLY DISPERSED 2 INTO THE SINUS CAVITY DURING TRANSCRESTAL MAXILLARY SINUS FLOOR ELEVATION: TWO-CASE REPORT
SINUSITIS, 2021, 5, 132-140

LOPEZ MA, PASSARELLI PC, RELLA E, NETTI A, LOPEZ A, CASALE M, D'ADDONA A
ALVEOLAR RIDGE AUGMENTATION WITH THE BONE INTO BONE TECHNIQUE: A HISTOLOGICAL AND HISTOMORPHOMETRIC ANALYSIS
J OSSEOINTEGR, SEPTEMBER 2021; 13(3)

BELTRAMINI G, RUSSILLO A, BASERGA C, PELLATI A, PIVA A, CANDOTTO V
HETEROLOGOUS CORTICO-CANCELLEUS BONE MIX STIMULATED DENTAL PULP DERIVED STEM CELLS
J BIOL REGUL HOMEOST AGENTS, JAN-FEB 2020;34(1 SUPPL. 2): 1-5. DENTAL SUPPLEMENT

DUMITRESCU CR, NEACSU IA, SURDU VA, NICOARA AI, IORDACHE F, TRUSCA R, CIOCAN LT, FICAI A, ANDRONESCU E
NANO-HYDROXYAPATITE VS. XENOGRAFTS: SYNTHESIS, CHARACTERIZATION, AND IN VITRO BEHAVIOR
NANOMATERIALS (BASEL), 2021 SEP 2;11(9):2289

FOTI V, SAVIO D, ROSSI R
ONE-TIME CORTICAL LAMINA: A NEW TECHNIQUE FOR HORIZONTAL RIDGE AUGMENTATION. A CASE SERIES
J BIOMED ENG MEDICAL IMAG, DECEMBER, 25, 2021, VOL. 8, NO. 6

ATIEH MA, ALSABEHA NH, PAYNE AG, ALI S, FAGGION CMJ, ESPOSITO M
INTERVENTIONS FOR REPLACING MISSING TEETH: ALVEOLAR RIDGE PRESERVATION TECHNIQUES FOR DENTAL IMPLANT SITE DEVELOPMENT
COCHRANE DATABASE SYST REV, 2021 APR 26;4(4):CD010176

GRASSI A, BERNARDELLO F, CAVANI F, PALUMBO C, SPINATO S
THREE-PUNCH ALVEOLAR RIDGE RECONSTRUCTION TECHNIQUE: A NOVEL FLAPLESS APPROACH IN EIGHT CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, NOV-DEC 2021;41(6):875-884

TERHEYDEN H, MEIJER GJ, RAGHOEBAR GM
VERTICAL BONE AUGMENTATION AND REGULAR IMPLANTS VERSUS SHORT IMPLANTS IN THE VERTICALLY DEFICIENT POSTERIOR MANDIBLE: A SYSTEMATIC REVIEW AND META-ANALYSIS OF RANDOMIZED STUDIES
INT J ORAL MAXILLOFAC SURG, 2021 SEP;50(9):1249-1258

HERNÁNDEZ ALFARO F, GONZÁLEZ SOTO J, GIRALT HERNANDO M, MASÍA GRIDILLA J, RAFFAINI M, VALLS ONTANÓN A
INTERPOSITIONAL COLLAGENATED CANCELLOUS BONE BLOCKS FOR NASAL DORSUM AUGMENTATION: A NEW TECHNIQUE FOR NASOMAXILLARY HYPOPLASIA TREATMENT
J PLAST RECONSTR AESTHET SURG, 2021 JAN;74(1):223-243

LYSOKON YY, LUCHYNSKYI MA, ISKIV MO
MORPHOLOGICAL FEATURES OF BONE TISSUE REGENERATES IN EXPERIMENTAL ANIMALS WHEN USING OSTEOTROPIC DRUGS FOR RECOVERY BONE DEFECTS
J EDUCATION HEALTH SPORT, 3 OCTOBER 2021, T. 11, NR 10, S. 11-17

AZAB E, YOUSSEF AR
BIOCOMPATIBILITY EVALUATION OF HUMAN AND PORCINE ACELLULAR DERMAL MATRIX ON HUMAN PRIMARY GINGIVAL FIBROBLASTS: IN VITRO COMPARATIVE STUDY
EUR J DENT, 2021 JUL;15(3):563-567

BEVILACQUA L, FONZAR A, OLIVIER S, DE BIASI M, VISINTIN M, ANGERAME D, MAGLIONE M
OUTCOME OF DIFFERENT SURGICAL APPROACHES IN THE TREATMENT OF CLASS II FURCATION DEFECTS IN MANDIBULAR MOLARS: A RANDOMIZED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, SEP/OCT 2020;40(5):693-701

BARGH S, SILINDIR-GUNAY M, OZER AY, PALASKA E, KARAARSLAN D, IDE S, SOLPAN D
PHYSICO-CHEMICAL EVALUATION OF GAMMA AND MICROWAVE IRRADIATED DENTAL GRAFTS
RADIATION PHYS CHEM, VOLUME 170, MAY 2020, 108627

BARGH S, SILINDIR-GUNAY M, OZER AY, COLAK S, KUTLU B, NOHUTCU R
THE EFFECTS OF GAMMA AND MICROWAVE STERILIZATION ON PERIODONTOLOGICAL GRAFTS
CHEM PHYS IMPACT, VOLUME 3, DECEMBER 2021, 100046

BARGH S, SILINDIR-GUNAY M, OZER AY, EKIZOGLU M, KART D, KUTLU B, NOHUTCU R
THE EFFECT OF GAMMA AND MICROWAVE RADIATION STERILIZATION ON PERIODONTOLOGICAL GRAFTS FOR MICROBIOLOGICAL EVALUATION
APMIS, 2021 OCT;129(10):587-597

CANELLAS JVDS, SOARES BN, RITTO FG, VETTORE MV, VIDIGAL JÚNIOR GM, FISCHER RG, MEDEIROS PJD
WHAT GRAFTING MATERIALS PRODUCE GREATER ALVEOLAR RIDGE PRESERVATION AFTER TOOTH EXTRACTION? A SYSTEMATIC REVIEW AND NETWORK METAANALYSIS
J CRANIOMAXILLOFAC SURG, 2021 NOV;49(11):1064-1071

TUMEDEI M, MANCINELLI R, DI FILIPPO ES, MARRONE M, IEZZI G, PIATTELLI A, FULLE S
OSTEOGENIC POTENTIAL OF HUMAN DENTAL PULP STEM CELLS CO-CULTURED WITH EQUINE BONE SUBSTITUTE COMBINED WITH MELATONIN
INT J PERIODONTICS RESTORATIVE DENT, JAN-FEB 2022;42(1):75-81

DHANDE S, ELASKARY A, MUGLIKAR S, SHEIKH S, HEGDE R, KALE R, JABADE B, SAYYED S
TREAT AND RESTORE COMPROMISED SOCKETS IMMEDIATELY WITH THE NOVEL VESTIBULAR SOCKET THERAPY: A REVOLUTION IN IMMEDIATE IMPLANT PLACEMENT
GLOBAL J MED RES: J DENTISTRY & OTOLARYNGOLOGY, VOLUME 21 ISSUE 3 VERSION 1.0 YEAR 2021

SALAH RM, ELBA GM, KAWANA KY
EFFECT OF COLLAGEN MEMBRANE ON HEALING OF ALVEOLAR BONE IN RATS WITH LIGATURE INDUCED PERIODONTITIS
ALEXANDRIA DENTAL J, AVAILABLE ONLINE FROM 25 FEBRUARY 2022

LYSOKON Y, BANDRIVSKY YL, LUCHYNSKYI MA
ANALYSIS OF THE RESULTS OF TREATMENT OF DESTRUCTIVE FORMS OF APICAL PERIODONTITIS WITH OSTEOTROPIC DRUGS IN A SHORT TERM
WIAD LEK, 2022;75(1 PT 2):228-231

CACCIANIGA G, REY G, BALDONI M, CACCIANIGA P, PORCARO G, BALDONI A, CERAULO S
LASER DECONTAMINATION AND LED PHOTOBIO-MODULATION PROMOTE BONE REGENERATION AND WOUND HEALING BY SECONDARY INTENTION, IN ALVEOLAR RIDGE PRESERVATION - CLINICAL AND RADIOGRAPHIC EVALUATION: A PILOT EXPERIENCE
PHOTOBIO-MODUL PHOTOMED LASER SURG, 2022 MAR 7

PERINI A, VIÑA ALMUNIA J, CARDA C, MARTÍN DE LLANO JJ, BOTTICELLI D, PEÑARROCHA DIAGO M
INFLUENCE OF THE USE OF A COLLAGEN MEMBRANE PLACED ON THE BONE WINDOW AFTER SINUS FLOOR AUGMENTATION - AN EXPERIMENTAL STUDY IN RABBITS
DENT J (BASEL), 2021 NOV 12;9(11):131

ELBARBARY A, REDA A, ELAZIZ AA
EVALUATION OF THE ADDITION OF INJECTABLE PLATELET RICH FIBRIN TO XENOGRAFT IN MANAGEMENT OF PERIODONTAL INTRAOSSEOUS DEFECTS. "RANDOMIZED CONTROLLED TRIAL"
AL-AZHAR DENTAL J, VOL. 9, NO. 2, APRIL (2022) - P. 321:330

GÜLER B, URAZ A, ÇETINER D
THE CHEMICAL EVALUATION OF DIFFERENT DENTAL GRAFT MATERIALS BY ENERGY DISPERSIVE X-RAY SPECTROMETRY TECHNIQUE
INT DENTAL RES, VOLUME 11, NUMBER 1, 2021

CAPELLA MONSONÍS H, ZEUGOLIS DI
DECELLULARIZED XENOGRAFTS IN REGENERATIVE MEDICINE: FROM PROCESSING TO CLINICAL APPLICATION
XENOTRANSPLANTATION, 2021 JUL;28(4):E12683

CORREIA F, POZZA DH, GOUVEIA S, CAMPOS FELINO A, FARIA ALMEIDA R
ADVANTAGES OF PORCINE XENOGRAFT OVER AUTOGRAFT IN SINUS LIFT: A RANDOMISED CLINICAL TRIAL
MATERIALS (BASEL), 2021 JUN 21;14(12):3439

PARUSHEVA S, METODIEV D, FAKIH H, BOSHPAKOVA T, GEORGIEV B
HISTOMORPHOLOGICAL STUDY OF REGENERATIVE PROCESSES IN POST-EXTRACTION ALVEOLAR SOCKETS AFTER APPLICATION OF A COMBINED PRESERVATION METHOD
C. R. ACAD. BULG. SCI., 74, NO 8, 2021

COMUZZI L, TUMEDEI M, PIATTELLI A, TARTAGLIA G, DEL FABBRO M
RADIOGRAPHIC ANALYSIS OF GRAFT DIMENSIONAL CHANGES IN TRANSCRESTAL MAXILLARY SINUS AUGMENTATION: A RETROSPECTIVE STUDY
MATERIALS (BASEL), 2022 APR 19;15(9):2964

COMUZZI L, TUMEDEI M, PIATTELLI A, TARTAGLIA G, DEL FABBRO M
RADIOGRAPHIC ANALYSIS OF GRAFT DIMENSIONAL CHANGES AFTER LATERAL MAXILLARY SINUS AUGMENTATION WITH HETEROLOGOUS MATERIALS AND SIMULTANEOUS IMPLANT PLACEMENT: A RETROSPECTIVE STUDY IN 18 PATIENTS
MATERIALS (BASEL), 2022 APR 22;15(9):3056

PISTILLI R, CANULLO L, PESCE P, PISTILLI V, CAPONIO C, SBRICOLI L
GUIDED IMPLANT SURGERY AND SINUS LIFT IN SEVERELY RESORBED MAXILLAE: A RETROSPECTIVE CLINICAL STUDY WITH UP TO 10 YEARS OF FOLLOW-UP
J DENT. 2022 JUN;121:104137

FELETTA L, BOTTICELLI D, APAZA ALCAYHUAMAN K, PEÑARROCHA DIAGO M, EZZEDDIN AYOUB M, ZARAGOZI ALONSO R, VIÑA ALMUNIA J
INFLUENCE OF THE USE OF AUTOGENOUS BONE PARTICLES TO CLOSE THE ACCESS WINDOW AFTER MAXILLARY SINUS FLOOR AUGMENTATION: A MICRO-COMPUTED TOMOGRAPHY AND POSITRON EMISSION TOMOGRAPHY STUDY IN RABBITS
ORAL MAXILLOFAC SURG (2022)

ESKAN MA
IMPLANT-ASSOCIATED GIANT CELL GRANULOMA: A CASE REPORT OF 4.8-YEAR FOLLOW-UP AND LITERATURE REVIEW
SELCUK DENTAL JOURNAL, 2022; 9

GURBUZ E, DURSUN E, VATANSEVER A, CAGLAYAN F
MICROCOMPUTED TOMOGRAPHIC ANALYSIS OF BONE MICROARCHITECTURE AFTER SINUS AUGMENTATION WITH HYALURONIC MATRIX: A CASE-CONTROL STUDY
ORAL MAXILLOFAC SURG. 2022 SEP;26(3):431-437

DI TINCO R, CONSOLO U, PISCIOTTA A, ORLANDI G, BERTANI G, NASI M, BERTACCHINI J, CARNEVALE G
CHARACTERIZATION OF DENTAL PULP STEM CELLS RESPONSE TO BONE SUBSTITUTES BIOMATERIALS IN DENTISTRY
POLYMERS, 2022 MAY 30;14(11):2223

MORIMOTO A, KOBAYASHI N, FERRI M, IEZZI G, PIATTELLI A, FORTICH MESA N, BOTTICELLI D
INFLUENCE ON IMPLANT BONE HEALING OF A COLLAGEN MEMBRANE PLACED SUBJACENT THE SINUS MUCOSA - A RANDOMIZED CLINICAL TRIAL ON SINUS FLOOR ELEVATION
DENT J (BASEL), 2022 JUN 8;10(6):105

TUMEDEI M, MIJIRITSKY E, MOURÃO CF, PIATTELLI A, DEGIDI M, MANGANO C, IEZZI G
HISTOLOGICAL AND BIOLOGICAL RESPONSE TO DIFFERENT TYPES OF BIOMATERIALS: A NARRATIVE SINGLE RESEARCH CENTER EXPERIENCE OVER THREE DECADES
INT J ENVIRON RES PUBLIC HEALTH, 2022, 19, 7942

TETÈ S, ZIZZARI VL, VINCI V, ZARA S, DI TORE U, MANICA M, CATALDI A, MORTELLARO C, PIATTELLI A, GHERLONE E
EQUINE AND PORCINE BONE SUBSTITUTES IN MAXILLARY SINUS AUGMENTATION: A HISTOLOGICAL AND IMMUNOHISTOCHEMICAL ANALYSIS OF VEGF EXPRESSION
J CRANIOFAC SURG, 2014 MAY;25(3):835-9

GHALLAB NA, ELASKARY A, ELSABAGH H, TOUKHY AE, ABDELRAHMAN H, EL-KIMARY G
A NOVEL ATRAUMATIC EXTRACTION TECHNIQUE USING VESTIBULAR SOCKET THERAPY FOR IMMEDIATE IMPLANT PLACEMENT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
ORAL MAXILLOFAC SURG. 2022 JUN 20

ELASKARY A, ELFANA A, MEABED M, ABD-ELWAHAB RADY I, AKRAM M, FAWZY EL-SAYED K
IMMEDIATE IMPLANT PLACEMENT UTILIZING VESTIBULAR SOCKET THERAPY VERSUS EARLY IMPLANT PLACEMENT WITH CONTOUR AUGMENTATION FOR REHABILITATION OF COMPROMISED EXTRACTION SOCKETS IN THE ESTHETIC ZONE: A RANDOMIZED CONTROLLED CLINICAL TRIAL
CLIN IMPLANT DENT RELAT RES. 2022 OCT;24(5):559-568

LOMBARDI T, LAMAZZA L, BERNARDELLO F, ZIĘTEK G, STACCHI C, TROIANO G
CLINICAL AND RADIOGRAPHIC OUTCOMES FOLLOWING TRANSCRESTAL MAXILLARY SINUS FLOOR ELEVATION WITH INJECTABLE XENOGENOUS BONE SUBSTITUTE IN GEL FORM: A PROSPECTIVE MULTICENTER STUDY
INT J IMPLANT DENT, 2022 JUL 22;8(1):32

ELASKARY A, ABDELRAHMAN H, ELSABAGH HH, EL-KIMARY GI
DOES GRAFTING THE JUMPING GAP IN IMMEDIATELY PLACED ANTERIOR IMPLANTS USING VESTIBULAR SOCKET THERAPY INFLUENCE THE LABIAL BONE THICKNESS?
J ORAL MAXILLOFAC SURG, 2022 AUG;80(8):1398-1407

OSORIO R, ASADY S, TOLEDANO OSORIO M, TOLEDANO M, BUENO JM, MARTÍNEZ OJEDA RM, OSORIO E
BIOMIMETIC REMINERALIZATION OF AN EXTRACELLULAR MATRIX COLLAGEN MEMBRANE FOR BONE REGENERATION
POLYMERS (BASEL), 2022 AUG 11;14(16):3274

ROMASCO T, TUMEDEI M, INCHINGOLO F, PIGNATELLI P, MONTESANI L, IEZZI G, PETRINI M, PIATTELLI A, DI PIETRO N
A NARRATIVE REVIEW ON THE EFFECTIVENESS OF BONE REGENERATION PROCEDURES WITH OSTEOBIOL® COLLAGENATED PORCINE GRAFTS: THE TRANSLATIONAL RESEARCH EXPERIENCE OVER 20 YEARS
J FUNCT BIOMATER, 2022 AUG 18;13(3):121

SIMONČIĆ B, DIMITRIJEVIĆ M, JULOSKI J, VULIČEVIĆ ZR
THE USE OF ADIPOSE-DERIVED STEM CELLS, PLATELET-RICH AND PLATELET-POOR PLASMA IN THE MAXILLARY CYST TREATMENT
SERBIAN ARCHIVES OF MEDICINE

IMAI H, IEZZI G, PIATTELLI A, FERRI M, APAZA ALCAYHUAMAN KA, BOTTICELLI D
INFLUENCE OF THE DIMENSIONS OF THE ANTROSTOMY ON OSSEOINTEGRATION OF MINI-IMPLANTS PLACED IN THE GRAFTED REGION AFTER SINUS FLOOR ELEVATION: A RANDOMIZED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS. 2020 MAY/JUN;35(3):591-598

MASUDA K, SILVA ER, BOTTICELLI D, ALCCAYHUAMAN KA, XAVIER SP
ANTROSTOMY PREPARATION FOR MAXILLARY SINUS FLOOR AUGMENTATION USING DRILLS OR A SONIC INSTRUMENT: A MICROCOMPUTED TOMOGRAPHY AND HISTOMORPHOMETRIC STUDY IN RABBITS
INT J ORAL MAXILLOFAC IMPLANTS. 2019 JULY/AUGUST;34(4):819-827

KANAYAMA M, BOTTICELLI D, ALCCAYHUAMAN K, YONEZAWA D, SILVA ER, XAVIER SP
THE IMPACT ON HEALING OF BIOACTIVATION WITH ARGON PLASMA OF A XENOGENIC GRAFT WITH ADEQUATE FIXATION BUT POOR ADAPTATION TO THE RECIPIENT SITE: AN EXPERIMENTAL STUDY IN RABBIT
INT J ORAL MAXILLOFAC IMPLANTS. 2021 JUL-AUG;36(4):703-714

LOPEZ MA, PASSARELLI PC, NETTI A, D'ADDONA A, CARINCI F, WYCHOWAŃSKI P, CECCHETTI F
PORCINE CORTICAL BONE LAMINA AS A PREDICTABLE TECHNIQUE FOR GUIDED BONE REGENERATION: HISTOMORPHOMETRIC AND RADIOGRAPHIC EVALUATION
APPL. SCI. 2022, 12, 10285

MALTAGLIATI A, ANGIERO F, FERRANTE F, BLASI S, OTTONELLO A
ALVEOLAR RIDGE EXPANSION BY IMPLANTS WITH OSTEODISTRACTIVE FUNCTION: A CLINICAL REPORT
DENT ORAL CRANIOFAC RES. 2016, VOLUME 2, ISSUE 6

FELICE P, LIZIO G, BARAUSSE C, ROCCOLI L, BONIFAZI L, PISTILLI R, SIMION M, PELLEGRINO G
REVERSE GUIDED BONE REGENERATION (R-GBR) DIGITALWORKFLOW FOR ATROPHIC JAWS REHABILITATION
APPL. SCI. 2022, 12,9947

ALQUTUB M, MUKHTAR A, ALALI Y, VOHRA F, ABDULJABBAR T
OSTEOGENIC DIFFERENTIATION OF PERIODONTAL LIGAMENT STEM CELLS SEEDED ON EQUINE-DERIVED XENOGRAFT IN OSTEOGENIC GROWTH MEDIA
MEDICINA (KAUNAS) 2022 OCT 25;58(11):1518

FOTI V, ROSSI R
FIBRINOGEN-INDUCED REGENERATION SEALING TECHNIQUE (F.I.R.S.T.): AN IMPROVEMENT AND MODIFICATION OF TRADITIONAL GBR: A REPORT OF TWO CASES
MOD RES DENT. 2020, 5(2)

WANIS RW, HOSNY MM, ELNAHASS H
CLINICAL EVALUATION OF THE BUCCAL ASPECT AROUND IMMEDIATE IMPLANT USING DUAL ZONE THERAPEUTIC CONCEPT VERSUS BUCCAL GAP FILL TO BONE LEVEL: A RANDOMIZED CONTROLLED CLINICAL TRIAL
CLIN IMPLANT DENT RELAT RES, 2022 JUN;24(3):307-319

KÖRNER G, BÄUMER-KÖNIG A, BARZ-POPP F, TRIPODAKIS A
ESTHETIC IMPLANT SOLUTIONS IN THE PERIODONTALLY COMPROMISED ANTERIOR REGION. EXPANDING THE INDICATIONS FOR IMMEDIATE IMPLANT PLACEMENT WITH THE SOCKET REBUILDING TECHNIQUE (SRT)
INT J ESTHET DENT. 2022 NOV 25;17(4):408-422

MIYAUCHI Y, IZUTANI T, TERANISHI Y, IIDA T, NAKAJIMA Y, XAVIER SP, BABA S
HEALING PATTERNS PATTERNS OF NON-COLLAGENATED BOVINE AND COLLAGENATED PORCINE XENOGRAFTS USED FOR SINUS FLOOR ELEVATION: A HISTOLOGICAL STUDY IN RABBITS
J. FUNCT. BIOMATER. 2022 DEC 5;13(4):276

HIROTA A, IEZZI G, PIATTELLI A, FERRI M, TANAKA K, APAZA ALCAYHUAMAN KA, BOTTICELLI D
INFLUENCE OF THE POSITION OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ON THE HEALING OF MINI-IMPLANTS: A RANDOMIZED CLINICAL TRIAL
ORAL MAXILLOFAC SURG. 2020 SEP;24(3):299-308

GRASSI A, MEMÈ L, STRAPPA EM, MARTINI E, BAMBINI F
MODIFIED PERIOSTEAL INHIBITION (MPI) TECHNIQUE FOR EXTRACTION SOCKETS: A CASE SERIES REPORT.
APPL. SCI. 2022, 12, 12292

PELÁEZ-CRUZ P, LÓPEZ JORNET P, TATULLO M, PONS-FUSTER LÓPEZ E
EPIGALLOCATECHIN-3-GALLATE IMPROVES THE BIOCOMPATIBILITY OF BONE SUBSTITUTES IN DENTAL PULP STEM CELLS
ANN ANAT. 2023 FEB;246:152045

NAKAJIMA Y, BOTTICELLI D, DE ROSSI EF, FERREIRA BALAN V, PIRES GODOY E, SILVA ER, XAVIER SP
SCHNEIDERIAN MEMBRANE COLLATERAL DAMAGE CAUSED BY COLLAGENATED AND NON-COLLAGENATED XENOGRAFTS: A HISTOLOGICAL STUDY IN RABBITS
DENT J (BASEL). 2023 JAN 26;11(2):31

CORREIA F, GOUVEIA SA, POZZA DH, FELINO AC, FARIA-ALMEIDA R
A RANDOMIZED CLINICAL TRIAL COMPARING IMPLANTS PLACED IN TWO DIFFERENT BIOMATERIALS USED FOR MAXILLARY SINUS AUGMENTATION
MATERIALS (BASEL). 2023 JAN 31;16(3):1220

RADUNOVIC M, PAVIC A, IVANOVIC V, MILIVOJEVIC M, RADOVIC I, DI CARLO R, PILATO S, FONTANA A, PIATTELLI A, PETROVIC S
BIOCOMPATIBILITY AND ANTIBIOFILM ACTIVITY OF GRAPHENE-OXIDE FUNCTIONALIZED TITANIUM DISCS AND COLLAGEN MEMBRANES
DENT MATER. 2022 JUL;38(7):1117-1127

COMUZZI L, TUMEDEI M, PETRINI M, ROMASCO T, LORUSSO F, DE ANGELIS F, PIATTELLI A, TATULLO M, DI PIETRO N
CLINICAL AND RADIOLOGICAL EVALUATION OF A SELF-CONDENSING BONE IMPLANT IN ONE-STAGE SINUS AUGMENTATION: A 3-YEAR FOLLOW-UP RETROSPECTIVE STUDY
INT. J. ENVIRON. RES. PUBLIC HEALTH JAN 31;20(3):2583

PIERFELICE TV, D'AMICO E, D'ERCOLE S, LEPORE S, PIATTELLI A, BARONE A, IEZZI G, PETRINI M
FUNCTIONALIZATION OF A CORTICAL MEMBRANE WITH A PHOTODYNAMIC PROTOCOL
J. FUNCT. BIOMATER. 2023, 14, 133

SAPOZNIKOV L, HAIM D, ZAVAN B, SCORTECCI G, HUMPHREY MF
A NOVEL PORCINE DENTIN-DERIVED BONE GRAFT MATERIAL PROVIDES EFFECTIVE SITE STABILITY FOR IMPLANT PLACEMENT AFTER TOOTH EXTRACTION: A RANDOMIZED CONTROLLED CLINICAL TRIAL
CLIN ORAL INVESTIG. 2023 FEB 24

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