**IK**

**TAPERED SHAPE  
INTERNAL HEXAGON**

# company profile

**RESISTA**<sup>®</sup>

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**RESISTA**, trademark owned by Ing. Carlo Alberto Issoglio & C. S.r.l. , identifies the full range of products designed by the company to meet the needs of dentists and dental technicians.

**RESISTA is known all over the world since 1946, and synonymous of guarantee, balance and continuity.**

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The company has one of the most ip-to-date manufacturing site for the medical devices, located in Omegna, a beautiful town on Orta Lake, Italy.

With commitment, merit and pride, the technical and scientific staff of the R&D department is daily involved to promote innovation and the improvement of our products and processes.

Over the years, the **RESISTA** brand has also become synonymous of fast and qualitative response to market changes and the needs of modern dentistry.

Quality, Price and Service summarize in three simple words the daily perception of the users of **RESISTA** branded products.

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**RESEARCH & DEVELOPMENT** - The R&D Resista Group's program stems from indications from the clinical world combined with our experience in the field of implantable medical devices.

Due to the use of 3D modeling programs and sophisticated computer systems, we simulate the final geometries and the design, developing the rapid prototyping phases of the projects.

**INNOVATIONS** - Resista Group is structured to meet all the needs of the metal-implant with an advanced conception technology, prototyping and final product engineering to arrive at innovative solutions. Validation is carried out in collaboration with accredited Research Institutes following the International Standards Requirements.

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**PROJECT AND PROCESS QUALITY** - ICIM Spa certified **Ing. C. A. Issoglio & C. S.r.l.** in accordance with **UNI EN ISO 9001** and **UNI EN ISO 13485** in compliance with all current regulations relating to the products and services offered. Every employee covering different roles (technical, engineering, administrative, commercial, mechanic operators, etc.), follows guidelines and objectives for a single purpose: continuous product improvement.

**MAIN OBJECTIVE** - Customer satisfaction is our objective. The winning factor is the company's ability to quickly resolve any type of request. The technical support provided by Product Specialists, availability, efficiency and friendliness are the strengths of our structure.

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**BIOMECHANICAL CONCEPTS** - The **Resista Group** method success is the result of several product design studies and the high technology applied in the production. The precision and reliability of the implant has a strong influence on the long term clinical success.


**The manufacturing process** is carried out by high qualified operators, skilled enough to obtain optimal and reproducible results with the use of electronically controlled CNC machinery. The machinery tolerances, especially in **prosthetic connections**, are our **gold standard** and they are made between **7 - 10 µm**.

**MATERIALS** - The implants are made in **Titanium ASTM Gr.4 Cold Worked** (ISO 5832/2), the prosthetic screws and the prosthetic components are made in **Titanium ASTM Gr.5 ELI** (ISO 5832/2).

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**SURFACE TREATMENTS** - The non-contaminant micro-subtraction process, **DAE** (Double Acid Etching), changes the implant micro roughness (**Ra, Rq**) and the surface texture maximizing the wettability and biomimetic properties. The new surface treatment and decontamination process are patented by **Nobil Bio Ricerche** improving the bio-chemical properties of the implants.


**QUALITY SYSTEM (QS)** - **UNI EN ISO 9001** and **UNI EN ISO 13485** in compliance of Medical Devices Directives. We are qualified in the design and in the production management of dental implants, dental prosthesis, intraligamental anesthesia syringes and abrasive discs.

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
**TITANIUM GR 5 PRE-ASSEMBLED MOUNTER** - One of the strong points of our **Resista** line is the Multifunction Moulder Device that is used in the **4** following different applications:

**1** Implant pick up and screwing - **2** Closed Tray Impression Transfer - **3** Definitive Abutment up to 13° angular correction - **4** Moulder for Computer Guided Surgery

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
**TRANSFER MOUNTER** - Transfer for **Closed Tray Technique** (Elastic tear impression technique).

It is a simple, fast and effective technique, suitable for all types of prosthetic projects. The Transfer Mounter can be repositioned in the impression only in one position, defined by the different sides and retentions arranged on the surface.


**ABUTMENT MOUNTER** - Definitive Prosthetic Abutment complete with definitive through screw.

The Titanium Gr 5 Abutment Mounter can be milled in height and angle up to 13°, using laboratory drills for Titanium.

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
**GCS MOUNTER** - Mount suitable for implants placement by **Computer Guided Surgery**.


The cylindrical diameter of the Moulder glides on the axis inside the Bushings for Computer Guided Surgery, whereas the 6 hexagonal alignment reference points (notches) define the final position.

**MOUNTER FUNCTION INSIDE THE PACKAGING** - Hangs the implant inside the packaging.





The **Titanium Bushing** contains the Moulder and the implant and guarantees the quality of the surface during the removal surgical step, avoiding the contamination due to the scratching.

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








## The Micro-Nano Rough DAE surface treatment speed up bone healing processes

1. Removes the manufacturing organic residuals along with Argon Plasma Cleaning
2. Smooths the edge and eliminates the micro-defects
3. Increases surface and wettability (hydrophilicity surface), improving the first fibrin bridges adhesion
4. Increases protein adhesion
5. Maximizes cellular adhesion through the Micro-Nano roughness suitable for actin filaments (phyllopods) anchoring
6. Changes the Titanium surface chemical characteristics improving the biocompatibility and increasing the cellular proliferation and vitality

## Surface topography and chemistry interact with the cellular differentiation processes


The implant surface uses the microtopography as a communication language with the host tissue cells. The Plasma Argon Cold Cleaning with a controlled-environment packaging at a strict protocol procedures, play a key-role in the control of adhered endotoxins (main immunological response player to implant surfaces).

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### VEGF/VEGF-R/RUNX2 Upregulation in Human Periodontal Ligament Stem Cells Seeded on Dual Acid Etched Titanium Disk

Francesca Diomede, Guya Diletta Marconi, Marcos F. X. B. Cavalcanti, Jacopo Pizzicannella, Sante Donato Pierdomenico, Luigia Fonticoli, Adriano Piattelli and Oriana Trubiani

Materials 2020, 13, 706; doi:10.3390/ma13030706

## SURFACE TREATMENT DISTRIBUTION ON IMPLANT'S BODY

All the **Resista's** implants have a surface treatment that is locally differentiated. The logic of surface treatment distribution is dictated by the endo-osseous implant position which requires different characteristics dependin on the muco-gingival component presence around.

1. **Machined Surface:** Brushed, rectify and protected during the roughening process
2. **Nano-rough Surface:** Speedy treatment, nano-rough with average Ra less than one micron
3. **Micro-rough Surface:** 3/5 $\mu$ : DAE Treatment, micro-rough with average Ra between 3/5 $\mu$
4. **Micro-rough Surface:** 4/7 $\mu$ : DAE Slow Treatment, micro-rough with average Ra between 4/7 $\mu$

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### The Bacterial Anti-Adhesive Activity of Double-Etched Titanium (DAE) as a Dental Implant Surface

Morena Petrini, Alessandra Giuliani, Emanuela Di Campi, Silvia Di Lodovico, Giovanna Iezzi, Adriano Piattelli and Simonetta D'Ercole<sup>1</sup>  
 International Journal of Molecular Sciences 2020, 21, 8315; doi:10.3390/ijms21218315

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## THE RESISTA'S IMPLANT LINE STRONG POINTS

- 1. Macro Geometry:** A vast range of implant shapes for different surgical needs (Cylinder, Tapered, Active, Short, Mini and Extra-Large)
- 2. Micro Geometry:** Last Generation Differentiated Surface (Full, Half and Machined) treated (DAE Micro and Nano rough), decontaminated in Argon Cold Plasma Reactor.
- 3. Multifunctional Mounter:** 4 components in a single article (Mounter, Transfer, Abutment and Mounter for Computer Guided Surgery) to maximize result and minimize costs.
- 4. Prosthetic Parts:** A very wide range of prosthetic variants with 7/10 micron machining tolerances on the connection, with emerging curvilinear profiles, flat and conical connections, platform switching, golden titanium screws and PVD TiN coatings (more aesthetic and less retention of dental plaque).
- 5. Surgical Instruments:** Surgical Drills **3-Tech**, maximum efficiency, minimum trauma, smaller friction coefficient, **PRO MSD** technology for counterclockwise bone densifying application long, lasting resistance, perfect visibility of laser marking and surgical millimeters stop.

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
 **Various bio-mechanical factors affecting heat generation during osteotomy preparation: A systematic review**

Chirag J Chauhan<sup>1</sup>, Darshana N Shah<sup>1</sup>, Foram B Sutaria<sup>1</sup>

Indian J Dent Res. Jan-Feb 2018;29(1):81-92. doi: 10.4103/ijdr.IJDR\_729\_16.

## LOGIC RELATION BETWEEN DRILL AND IMPLANT CORE

The rotary instruments design and shape is based on the logical coherence between the implants size, in the external dimensions, internal core and their penetration properties.

The **IC** line **Cylindrical Implants** have a tapered apical portion 3mm long.

The **IK** line **Conical Implants** have a tapered apical portion 5mm long.

The drills have respectively 2/3 apical sections, suitable for the implant tunnel preparation, with respect to the implant geometry, maintaining the best bone implant contact (BIC).

Each step drill is equipped with a sharp cutting edge to facilitate penetration, minimizing friction and overheating.

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
### ✓ Bone classification: clinical-histomorphometric comparison

Trisi P, Rao W. - Clin Oral Implants Res. 1999 Feb;10(1):1-7. doi: 10.1034/j.1600-0501.1999.100101.x.

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## THE CONDITIONING OF SURGICAL PROTOCOLS

The differentiated implant geometries in the macro/micro shape provide for different applications, both according to the surgical/prosthetic requirements (immediate loading, early loading or postponed loading) and according to the conditions of the bone/anatomical component of the patient who needs rehabilitation. However, surgical maneuvers and alternative protocols are often required to improve implant performance, according to the morphological characteristics of the residual bone structures.

The success of these maneuvers, that is often the result of a customized tools and protocols combination, represents the different response to the problem in relation to the personal operating skills and the presence of tools suitable for their application.

## MECHANICAL OSSEODENSIFICATION IN ANTI-CLOCKWISE ROTATION

The osseodensification using rotary instruments is a recent surgical technique for the implant site preparation that can be associated with different protocols; these protocols can be used in those particular anatomical conditions such as poor bone quality and insufficient vertical / horizontal dimensions.

This non-subtractive osteotomic bone-condensation approach produces an increase in peri-osteotomic bone density, saving the tissue and increasing the primary implant stability.

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


## RESISTA PRO MSD SURGICAL DRILLS - Modular Surgical Drilling

The Modular Surgical Drills with variable geometry are 3-section helicoidal drills with progressive radius rake, that can be used clockwise and anticlockwise.


The diameter of each drill changes according to a constant ratio of 0.4mm (30/34/38/42/46/50): this allows the operator to choose the drill according to the bone quality (over-preparation or under-preparation). All the drills are coated with a DLC (Diamond Like Carbon) synthetic diamond solid lubricant that maximizes performance in terms of mechanical strength and friction reduction.

The **anticlockwise rotation**, reversing the involved forces, generates 3 different effects on the implant tunnel perforation, which can revolutionize the known milling logic.


- 1) Anterior and lateral thrust of the bone removed by the tip, and of the liquids as well, blood and physiological water.
- 2) Ejection and "hammer effect" of the drill that produce an improvement in vertical control.
- 3) Reduction of cutting efficiency to protect sensitive anatomical parts.

# macro geometry

## IK tapered FULL TREATMENT

## IK tapered HT HALF TREATMENT

# macro geometry

## IK tapered ML FULL TREATMENT

●  
3.4

●  
3.8

●  
4.2

●  
4.6

●  
5.0

## IK tapered XL ML FULL TREATMENT

●  
6.0


●  
8.0

# INTERNAL HEXAGON UNIVERSAL CONNECTION

TAPERED  
SHAPE

IK


IMPIANTO  
CON MOUNTER  
WITH MOUNTER


# IK 34

**K-TAPERED**  
internal hexagon - full treatment

## TAPERED SHAPE + MOUNTER

CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 3407	Ø 3.4	7.0 mm	Ø 3.5 CPS	Ø 1.8
IK 3408	Ø 3.4	8.5 mm	Ø 3.5 CPS	Ø 1.8
IK 3410	Ø 3.4	10 mm	Ø 3.5 CPS	Ø 1.8
IK 3411	Ø 3.4	11.5 mm	Ø 3.5 CPS	Ø 1.8
IK 3413	Ø 3.4	13 mm	Ø 3.5 CPS	Ø 1.8
IK 3415	Ø 3.4	15 mm	Ø 3.5 CPS	Ø 1.8


# K-TAPERED


internal hexagon - full treatment

# IK 38

## TAPERED SHAPE + MOUNTER


CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 3808	Ø 3.8	8.5 mm	Ø 4.0	Ø 1.9
IK 3810	Ø 3.8	10 mm	Ø 4.0	Ø 1.9
IK 3811	Ø 3.8	11.5 mm	Ø 4.0	Ø 1.9
IK 3813	Ø 3.8	13 mm	Ø 4.0	Ø 1.9
IK 3815	Ø 3.8	15 mm	Ø 4.0	Ø 1.9


# IK 42

**K-TAPERED**  
internal hexagon - full treatment

## TAPERED SHAPE + MOUNTER

CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 4207	Ø 4.2	7.0 mm	Ø 4.0	Ø 2.3
IK 4208	Ø 4.2	8.5 mm	Ø 4.0	Ø 2.3
IK 4210	Ø 4.2	10 mm	Ø 4.0	Ø 2.3
IK 4211	Ø 4.2	11.5 mm	Ø 4.0	Ø 2.3
IK 4213	Ø 4.2	13 mm	Ø 4.0	Ø 2.3
IK 4215	Ø 4.2	15 mm	Ø 4.0	Ø 2.3


# K-TAPERED


internal hexagon - full treatment

# IK 46

## TAPERED SHAPE + MOUNTER


CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 4607	Ø 4.6	7.0 mm	Ø 4.0	Ø 2.5
IK 4608	Ø 4.6	8.5 mm	Ø 4.0	Ø 2.5
IK 4610	Ø 4.6	10 mm	Ø 4.0	Ø 2.5
IK 4611	Ø 4.6	11.5 mm	Ø 4.0	Ø 2.5
IK 4613	Ø 4.6	13 mm	Ø 4.0	Ø 2.5
IK 4615	Ø 4.6	15 mm	Ø 4.0	Ø 2.5


# IK 34 HT

**K-TAPERED HT**  
internal hexagon - half treatment

## TAPERED SHAPE + MOUNTER

CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 3408 HT	Ø 3.4	8.5 mm	Ø 3.5 CPS	Ø 1.8
IK 3410 HT	Ø 3.4	10 mm	Ø 3.5 CPS	Ø 1.8
IK 3411 HT	Ø 3.4	11.5 mm	Ø 3.5 CPS	Ø 1.8
IK 3413 HT	Ø 3.4	13 mm	Ø 3.5 CPS	Ø 1.8
IK 3415 HT	Ø 3.4	15 mm	Ø 3.5 CPS	Ø 1.8


# K-TAPERED HT


internal hexagon - half treatment

# IK 38 HT

## TAPERED SHAPE + MOUNTER


CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 3808 HT	Ø 3.8	8.5 mm	Ø 4.0	Ø 1.9
IK 3810 HT	Ø 3.8	10 mm	Ø 4.0	Ø 1.9
IK 3811 HT	Ø 3.8	11.5 mm	Ø 4.0	Ø 1.9
IK 3813 HT	Ø 3.8	13 mm	Ø 4.0	Ø 1.9
IK 3815 HT	Ø 3.8	15 mm	Ø 4.0	Ø 1.9


# IK 42 HT

**K-TAPERED HT**  
internal hexagon - half treatment

## TAPERED SHAPE + MOUNTER

CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 4208 HT	Ø 4.2	8.5 mm	Ø 4.0	Ø 2.3
IK 4210 HT	Ø 4.2	10 mm	Ø 4.0	Ø 2.3
IK 4211 HT	Ø 4.2	11.5 mm	Ø 4.0	Ø 2.3
IK 4213 HT	Ø 4.2	13 mm	Ø 4.0	Ø 2.3
IK 4215 HT	Ø 4.2	15 mm	Ø 4.0	Ø 2.3


# K-TAPERED HT


internal hexagon - half treatment

# IK 46 HT

## TAPERED SHAPE + MOUNTER

CODE	IMPLANT	LENGTH	PLATFORM	APEX
IK 4608 HT	Ø 4.6	8.5 mm	Ø 4.0	Ø 2.5
IK 4610 HT	Ø 4.6	10 mm	Ø 4.0	Ø 2.5
IK 4611 HT	Ø 4.6	11.5 mm	Ø 4.0	Ø 2.5
IK 4613 HT	Ø 4.6	13 mm	Ø 4.0	Ø 2.5
IK 4615 HT	Ø 4.6	15 mm	Ø 4.0	Ø 2.5


# INTERNAL HEXAGON

**WIDE PLATFORM**

**DOUBLE CONNECTION**

LARGE TAPERED  
SHAPE

IK


IMPIANTO  
CON MOUNTER  
WITH MOUNTER


# IK 50

**K-TAPERED**  
internal hexagon - full treatment

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













**IP 4HA2-W**      H 2 mm - **CPS** Standard Gingival Profile  
**IP 4HA4-W**      H 4 mm - **CPS** Standard Gingival Profile  
**IP 4HA6-W**      H 6 mm - **CPS** Standard Gingival Profile

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**IP 4HA2-E**      H 2 mm - Large Gingival Profile  
**IP 4HA4-E**      H 4 mm - Large Gingival Profile  
**IP 4HA6-E**      H 6 mm - Large Gingival Profile

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Temporary Glueing Abutments - Titanium Gr. 5

**IP 4TA-W** Temporary Glueing Abutment - Antirotation Hex

**IP 4TR-W** Temporary Glueing Abutment - Rotating

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**DM 1/2** Standard Hexagonal Driver Ø 1.27

## Pilastro Temporaneo da Saldatura - Titanio Gr. 5

Temporary Welding Abutments - Titanium Gr. 5

**IP 4TRS** Temporary Welding Abutment - Rotating

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**DM 1/2** Standard Hexagonal Driver Ø 1.27

## Moncone Titanio Diritto - Titanio Gr. 5 - PVD TiN

Straight Titanium Abutment - Titanium Gr. 5 - PVD TiN

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**DM 1/2** Standard Hexagonal Driver Ø 1.27

## Moncone Titanio Angolato - Titanio Gr. 5 - PVD TiN


Angled Titanium Abutment - Titanium Gr. 5 - PVD TiN

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
**Standard Link for Bonding Technique - Titanium Gr. 5**

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**Flat Digital Link for Bonding Technique - Titanium Gr. 5**

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**DM 1/2** Standard Hexagonal Driver Ø 1.27


**Transfer per impronta Pick Up - Titanio Gr. 5**


**Pick Up Impression Coping - Titanium Gr. 5**

**IP X40-W** H 12mm - Internal Hexagon CPS

**DM 1/2** Standard Hexagonal Driver Ø 1.27

# double connection friction fit 1.5°

32Ncm

## Moncone Titanio Dritto - Titanio Gr. 5 - PVD TiN

**Straight Titanium Abutment** - Titanium Gr. 5 - PVD TiN

- IP 4ST2-WF** 0° - Straight H2mm - **CPS** Friction Fit
- IP 4ST3-WF** 0° - Straight H3.5mm - **CPS** Friction Fit
- IP 5ST2-F** 0° - Straight H2mm (Large) - **FTF** Friction Fit

**IP VITG**  
**DM 1/2**

Vite di Ricambio Gold - Gold Spare Screw  
Standard Hexagonal Driver Ø 1.27

32Ncm

## Moncone in Titanio Angolato 15° - Titanio Gr. 5 - PVD TiN

**15° Angled Titanium Abutment** - Titanium Gr. 5 - PVD TiN

- IP 4AT152-WF** 15° Angled H2mm - **CPS** Friction Fit
- IP 4AT153-WF** 15° Angled H3mm - **CPS** Friction Fit
- IP 5ST152-F** 15° - Angled H2mm (Large) - **FTF** Friction Fit

**IP VITG**  
**DM 1/2**

Vite di Ricambio Gold - Gold Spare Screw  
Standard Hexagonal Driver Ø 1.27

32Ncm

## Moncone in Titanio Angolato 25° - Titanio Gr. 5 - PVD TiN

**25° Angled Titanium Abutment** - Titanium Gr. 5 - PVD TiN

- IP 4AT252-WF** 25° Angled H2mm - **CPS** Friction Fit
- IP 4AT253-WF** 25° Angled H3mm - **CPS** Friction Fit
- IP 5ST252-F** 25° - Angled H2mm (Large) - **FTF** Friction Fit

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
## 45° CONICAL SHOULDER CLOSURE

The 45° internal prosthetic support, already used in internal Zimmer type connections since the 90s, has an indisputable bibliography and moreover guarantees the distribution of vector forces that minimize the lateral bending moment.

This type of conical support, normally used in the distribution mechanics, increases the contact surface in relation to a plane surface and releases the forces on a curve of multi-directional vectors.

## FRICTION FIT 1.5° CONICAL HEX CONNECTION

The 1.5° conicity on the flat surfaces of the male / female hexagon causes a so close reduction of the mechanical interface that confers a direct binding (known as "cold welding") between the abutment and the implant.

This effect is fully activated at the end of the prosthetic tightening of the passing screw, that reaches 30 Ncm of torque and guarantees a bacterial seal through the reduction of the interface gap to zero.

In order to remove the abutment it will be necessary to use an alternative removal screw, that needs to be screwed in the original one place: its unique function is to push the abutment from the inside vertically and to remove it from its bound position.

**Transfer per Impronta Digitale - Titanio Gr. 5 + MS**  
**Scan Abutment - Titanium Gr. 5 + Micro Sandblasting**

**IP 4SA2 DS**      Antirotation Hex - **Narrow** Gingival Profile

**DM 1/2**              Standard Hexagonal Driver Ø 1.27

**Analogo CAD da Laboratorio con Vite - Titanio Gr. 5**  
**Laboratory CAD Analog with Screw - Titanium Gr. 5**

**IP 4AL2 DS**      Antirotation Hex - Ø 4 Universal Shape

**DM 1/2**              Standard Hexagonal Driver Ø 1.27

**Link Base per Incollaggio - Titanio Gr. 5**  
**Bonding Technique Base Link - Titanium Gr. 5**

**IP 4TB0 DS-W**      Gingival Profile **H 0.3mm** - Antirotation Hex  
**IP 4TBR0 DS-W**      Gingival Profile **H 0.3mm** - Rotating

**IP VITG**              Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2**              Standard Hexagonal Driver Ø 1.27

Adjustable in 3 heights


**Link Base per Incollaggio - Titanio Gr. 5**  
**Bonding Technique Base Link - Titanium Gr. 5**

**IP 4TB1 DS-W**      Gingival Profile **H 1.5mm** - Antirotation Hex  
**IP 4TBR1 DS-W**      Gingival Profile **H 1.5mm** - Rotating

**IP VITG**              Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2**              Standard Hexagonal Driver Ø 1.27

Adjustable in 3 heights

**Link Base per Incollaggio - Titanio Gr. 5**  
**Bonding Technique Base Link - Titanium Gr. 5**


**IP 4TB2 DS-W**      Gingival Profile **H 2.5mm** - Antirotation Hex  
**IP 4TBR2 DS-W**      Gingival Profile **H 2.5mm** - Rotating

**IP VITG**              Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2**              Standard Hexagonal Driver Ø 1.27

Adjustable in 3 heights

# toronto universal prosthetics

## Pilastro Toronto Dritto - Titanio Gr. 5

Toronto Straight Abutment - Titanium Gr. 5

- IP MS1-W H 1 mm - CPS Standard Gingival Profile
- IP MS2-W H 2 mm - CPS Standard Gingival Profile
- IP MS3-W H 3 mm - CPS Standard Gingival Profile

**DCDM** Toronto Straight Abutment Driver

## Pilastro Toronto Angolato 17° - Titanio Gr. 5

17° Toronto Angled Abutment - Titanium Gr. 5

- IP MA172-W 17° - H 2 mm - CPS Standard Gingival Profile
- IP MA173-W 17° - H 3 mm - CPS Standard Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw  
**DMX 1/2** Standard Torx Driver

## Pilastro Toronto Angolato 30° - Titanio Gr. 5

30° Toronto Angled Abutment - Titanium Gr. 5

- IP MA303-W 30° - H 3 mm - CPS Standard Gingival Profile
- IP MA304-W 30° - H 4 mm - CPS Standard Gingival Profile
- IP MA305-W 30° - H 5 mm - CPS Standard Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw  
**DMX 1/2** Standard Torx Driver

**Pilastro Temporaneo Toronto - Titanio Gr. 5**  
**Toronto Temporary Abutment - Titanium Gr. 5**

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**Toronto Spare Screws** - Titanium Gr. 5

**EP MV**                      Prosthetic Spare Screw  
**EP MVL**                    Impression Spare Screw

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
**CPS Cover Screw** - Titanium Gr. 4

**IP 4CS1**                    CPS Cover Screw - 1.5mm  
**IP 4CS3**                    CPS Cover Screw - 3.0mm  
**IP 4CS5**                    CPS Cover Screw - 5.0mm

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







**Healing Abutment** - Titanium Gr. 5 - Anodized

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**Laboratory Analog** - Titanium Gr. 5 - Anodized

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**IP 5CR** Rotating - Large Gingival Profile

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32Ncm 

**Monconi in Titanio** - Titanio Gr. 5 - PVD TiN

**Titanium Abutments** - Titanium Gr. 5 - PVD TiN

**IP 5ST3** 0° - Straight - Large Gingival Profile

**IP 5AT15** 15° - Angled - Large Gingival Profile

**IP 5AT25** 25° - Angled - Large Gingival Profile

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**Ucla Gold Abutment** - Gold/Platinum/Iridium Base

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## Pilastrini Temporanei - Titanio Gr. 5

Temporary Abutments - Titanium Gr. 5

**IP 4TA-W** Temporary Glueing Abutment - Antirotation Hex

**IP 4TR-W** Temporary Glueing Abutment - Rotating

**IP 4TRS** Temporary Welding Abutment - Rotating

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Aesthetic Temporary Abutment - Medical PEEK

**IP 5PA** Antirotation Hex - Large Gingival Profile

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## Moncone Estetico - Zirconia + Link in Titanio Gr. 5

Aesthetic Abutment - Zirconium + Titanium Gr. 5 Link

**IP 5SZ** 0° - Straight - Large Gingival Profile

**IP 5AZ15** 15° - Angled - Large Gingival Profile

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32Ncm

## Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

**IP X50** H 12mm - Internal Hexagon

**DM 1/2** Standard Hexagonal Driver

25Ncm

# wide prosthetics

## Pilastro a Sfera per Overdenture - Titanio Gr. 5 - PVD TiN

Overdenture Ball Abutment - Titanium Gr. 5 - PVD TiN

32Ncm 

- IP 4BA1-W** H 1 mm - **CPS** Standard Gingival Profile - NORMO
- IP 4BA2-W** H 2 mm - **CPS** Standard Gingival Profile - NORMO
- IP 4BA3-W** H 3 mm - **CPS** Standard Gingival Profile - NORMO
- IP 4BA4-W** H 4 mm - **CPS** Standard Gingival Profile - NORMO
- IP 4BA6-W** H 6 mm - **CPS** Standard Gingival Profile - NORMO

**DM 1/2** [Standard Hexagonal Driver](#)

## Kit Pilastro Equator per Overdenture - Titanio Gr. 5 - PVD TiN

Overdenture Equator Abutment Kit - Titanium Gr. 5 - PVD TiN

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Flat Digital Link for Bonding Technique - Titanium Gr. 5

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
Laboratory Digital Analog - Titanium Gr. 5 - Anodized

**IP 5AL2DS** H 12 mm - Universal Shape


# toronto universal prosthetics

## Pilastro Toronto Diritto - Titanio Gr. 5

Toronto Straight Abutment - Titanium Gr. 5

- IP MS1-W H 1 mm - Standard Gingival Profile
- IP MS2-W H 2 mm - Standard Gingival Profile
- IP MS3-W H 3 mm - Standard Gingival Profile

**DCDM** Toronto Straight Abutment Driver

## Pilastro Toronto Angolato 17° - Titanio Gr. 5

17° Toronto Angled Abutment - Titanium Gr. 5

- IP MA172-W 17° - H 2 mm - Standard Gingival Profile
- IP MA173-W 17° - H 3 mm - Standard Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw  
**DMX 1/2** Standard Torx Driver

## Pilastro Toronto Angolato 30° - Titanio Gr. 5

30° Toronto Angled Abutment - Titanium Gr. 5

- IP MA303-W 30° - H 3 mm - Standard Gingival Profile
- IP MA304-W 30° - H 4 mm - Standard Gingival Profile
- IP MA305-W 30° - H 5 mm - Standard Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw  
**DMX 1/2** Standard Torx Driver

## Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

- EP MT Rotating - Bonding Technique - Incollaggio Passivo
- EP MTS Rotating - Welding Technique - Saldatura Intraorale

EP MV Vite di Ricambio Toronto - Toronto Spare Screw  
**DMX 1/2** Standard Torx Driver

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

Toronto Healing Abutment and Impression Coping




- EP MI Rotating - Standard Gingival Profile
- EP MHA H 4 mm - Standard Gingival Profile
- EP MHA6 H 6 mm - Large Gingival Profile

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

Angled 17°/30° MUA Abutment Screw Ø 1.8mm

**EP MV**

Toronto Prosthetic Spare Screw - Ø 1.4mm

**EP MVL**

Toronto Impression Spare Screw - Ø 1.4mm

## RICAMBI MOUNTERS IMPIANTI - IMPLANT MOUNTERS SPARE PARTS

**MTA Mounter Transfer Abutment - Titanio Gr. 5**

**MTA Mounter Transfer Abutment - Titanium Gr. 5**

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







# guided surgical drills

## FG C - CORTICAL DRILLS

- FGR C35 S** Corticall Drill - Ø 3.5 mm - ● Small
- FGR C40 S** Corticall Drill - Ø 4.0 mm - ● Small
- FGR C35 L** Corticall Drill - Ø 3.5 mm - ● Large
- FGR C40 L** Corticall Drill - Ø 4.0 mm - ● Large
- FGR C50 L** Corticall Drill - Ø 5.0 mm - ● Large

## FG 20 S - FIRST DRILLS ● Small Driver

- FGR 2008 S** First Drill - Ø 2.0 mm - Long 8,5 mm
- FGR 2010 S** First Drill - Ø 2.0 mm - Long 10 mm
- FGR 2011 S** First Drill - Ø 2.0 mm - Long 11,5 mm
- FGR 2013 S** First Drill - Ø 2.0 mm - Long 13 mm
- FGR 2015 S** First Drill - Ø 2.0 mm - Long 15 mm

## FG 20 L - FIRST DRILLS ● Large Driver

- FG 2008 L** First Drill - Ø 2.0 mm - Long 8,5 mm
- FGR 2010 L** First Drill - Ø 2.0 mm - Long 10 mm
- FGR 2011 L** First Drill - Ø 2.0 mm - Long 11,5 mm
- FGR 2013 L** First Drill - Ø 2.0 mm - Long 13 mm
- FGR 2015 L** First Drill - Ø 2.0 mm - Long 15 mm

# guided surgical drills

## FG 34 L - 3.4 IMPLANT DRILLS ● Small Driver

- FGR 3408 S** 3.4 Implant Drill - Ø 2.9 mm - Long 8,5 mm
- FGR 3410 S** 3.4 Implant Drill - Ø 2.9 mm - Long 10 mm
- FGR 3411 S** 3.4 Implant Drill - Ø 2.9 mm - Long 11,5 mm
- FGR 3413 S** 3.4 Implant Drill - Ø 2.9 mm - Long 13 mm
- FGR 3415 S** 3.4 Implant Drill - Ø 2.9 mm - Long 15 mm

## FG 34 L - 3.4 IMPLANT DRILLS ● Large Driver

- FGR 3408 L** 3.4 Implant Drill - Ø 2.9 mm - Long 8,5 mm
- FGR 3410 L** 3.4 Implant Drill - Ø 2.9 mm - Long 10 mm
- FGR 3411 L** 3.4 Implant Drill - Ø 2.9 mm - Long 11,5 mm
- FGR 3413 L** 3.4 Implant Drill - Ø 2.9 mm - Long 13 mm
- FGR 3415 L** 3.4 Implant Drill - Ø 2.9 mm - Long 15 mm

## FG 38 S - 3.8 IMPLANT DRILLS ● Small Driver

- FGR 3808 S** 3.8 Implant Drill - Ø 3.3 mm - Long 8,5 mm
- FGR 3810 S** 3.8 Implant Drill - Ø 3.3 mm - Long 10 mm
- FGR 3811 S** 3.8 Implant Drill - Ø 3.3 mm - Long 11,5 mm
- FGR 3813 S** 3.8 Implant Drill - Ø 3.3 mm - Long 13 mm
- FGR 3815 S** 3.8 Implant Drill - Ø 3.3 mm - Long 15 mm

## FG 38 L - 3.8 IMPLANT DRILLS ● Large Driver

- FGR 3808 L** 3.8 Implant Drill - Ø 3.3 mm - Long 8,5 mm
- FGR 3810 L** 3.8 Implant Drill - Ø 3.3 mm - Long 10 mm
- FGR 3811 L** 3.8 Implant Drill - Ø 3.3 mm - Long 11,5 mm
- FGR 3813 L** 3.8 Implant Drill - Ø 3.3 mm - Long 13 mm
- FGR 3815 L** 3.8 Implant Drill - Ø 3.3 mm - Long 15 mm

# guided surgical drills

## FG 42 S - 4.2 IMPLANT DRILLS ● Small Driver

- FGR 4208 S** 4.2 Implant Drill - Ø 3.6 mm - Long 8,5 mm
- FGR 4210 S** 4.2 Implant Drill - Ø 3.6 mm - Long 10 mm
- FGR 4211 S** 4.2 Implant Drill - Ø 3.6 mm - Long 11,5 mm
- FGR 4213 S** 4.2 Implant Drill - Ø 3.6 mm - Long 13 mm
- FGR 4215 S** 4.2 Implant Drill - Ø 3.6 mm - Long 15 mm

## FG 42 L - 4.2 IMPLANT DRILLS ● Large Driver

- FGR 4208 L** 4.2 Implant Drill - Ø 3.6 mm - Long 8,5 mm
- FGR 4210 L** 4.2 Implant Drill - Ø 3.6 mm - Long 10 mm
- FGR 4211 L** 4.2 Implant Drill - Ø 3.6 mm - Long 11,5 mm
- FGR 4213 L** 4.2 Implant Drill - Ø 3.6 mm - Long 13 mm
- FGR 4215 L** 4.2 Implant Drill - Ø 3.6 mm - Long 15 mm

## FG 46 L - 4.6 IMPLANT DRILLS ● Large Driver

- FGR 4608 L** 4.6 Implant Drill - Ø 3.9 mm - Long 8,5 mm
- FGR 4610 L** 4.6 Implant Drill - Ø 3.9 mm - Long 10 mm
- FGR 4611 L** 4.6 Implant Drill - Ø 3.9 mm - Long 11,5 mm
- FGR 4613 L** 4.6 Implant Drill - Ø 3.9 mm - Long 13 mm
- FGR 4615 L** 4.6 Implant Drill - Ø 3.9 mm - Long 15 mm

## FG 50 L - 5.0 IMPLANT DRILLS ● Large Driver

- FGR 5008 L** 5.0 Implant Drill - Ø 4.3 mm - Long 8,5 mm
- FGR 5010 L** 5.0 Implant Drill - Ø 4.3 mm - Long 10 mm
- FGR 5011 L** 5.0 Implant Drill - Ø 4.3 mm - Long 11,5 mm
- FGR 5013 L** 5.0 Implant Drill - Ø 4.3 mm - Long 13 mm

# guided surgical drills

## GUIDED MOUNTERS

- IP4MG 4S** Mounter Guidata - Internal Hex ● Small Drive
- IP4MG 4L** Mounter Guidata - Internal Hex ● Large Drive
  
- SP3MG 3S** Mounter Guidata - External Hex ● Small Drive
- SP4MG 4S** Mounter Guidata - External Hex ● Small Drive
- SP4MG 4L** Mounter Guidata - External Hex ● Large Drive
  
- IP3MG 3S-P** Mounter Guidata - Platinum CC ● Small Drive
- IP3MG 3L-P** Mounter Guidata - Platinum CC ● Large Drive
- IP4MG 4L-P** Mounter Guidata - Platinum CC ● Large Drive

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


Kit OT Equator Internal Hexagon - Abutment + Box + Caps

<b>IP4EQ1</b>	Kit OT Equator h1 mm Ø4.0 CPS ●
<b>IP4EQ2</b>	Kit OT Equator h2 mm Ø4.0 CPS ●
<b>IP4EQ3</b>	Kit OT Equator h3 mm Ø4.0 CPS ●
<b>IP4EQ4</b>	Kit OT Equator h4 mm Ø4.0 CPS ●
<b>IP4EQ5</b>	Kit OT Equator h5 mm Ø4.0 CPS ●
<b>IP4EQ6</b>	Kit OT Equator h6 mm Ø4.0 CPS ●

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### Kit LOCATOR® Internal Hexagon - Abutment + Box + Caps

- IP4LC1** Kit LOCATOR h1 mm - Ø4.0 CPS ●  
**IP4LC2** Kit LOCATOR h2 mm - Ø4.0 CPS ●  
**IP4LC3** Kit LOCATOR h3 mm - Ø4.0 CPS ●  
**IP4LC4** Kit LOCATOR h4 mm - Ø4.0 CPS ●  
**IP4LC5** Kit LOCATOR h5 mm - Ø4.0 CPS ●  
**IP4LC6** Kit LOCATOR h6 mm - Ø4.0 CPS ●

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

**EXP** 4 Progressive diameters  
**KIT 005** Expander Surgical Kit

<b>A 1.3</b>	<b>A 1.6</b>	<b>A 1.9</b>	<b>A 2.4</b>
<b>B 2.1</b>	<b>B 2.8</b>	<b>B 3.5</b>	<b>B 4.3</b>
<b>C 2.4</b>	<b>C 3.1</b>	<b>C 3.8</b>	<b>C 4.6</b>

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
**PC 01R** Short - Hexagonal Implant Driver

**PC 02R** Medium - Hexagonal Implant Driver

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**AC 01** Short - Hexagonal Implant Driver

**AC 02** Medium - Hexagonal Implant Driver

# instruments and tools

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Manual Driver - Stainless Steel

**CRND** Cricchetto non Dinamometrico - No Torque Control Ratchet

CRD2-G

**CRD2** Cricchetto Dinamometrico 10/70 Ncm - Torque Control

**CRD2-G** Ghiera Cricchetto CRD2 con Attacco PC - CRD2 Ratchet Wheel with PC Connection

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








**KIT 007** Implant Extractor Kit - 2 Extractors + 1 Box + CRND

## Siringhe per Anestesia - Inox

### Anesthesia Syringes - Stainless Steel

PERIPEN

PERIPRESS

**22** PERIPRESS - Siringa Inox 1.8cc - Set 2 Portafiala

**32** PERIPRESS - Siringa Standard 1.8cc - Set 2 Portafiala

**118** PERIPEN - Siringa Inox 1.8cc - Set 2 Portafiala AM/EU

**Ing. C. A. Issoglio & C. S.r.l.**  
Via F.lli Di Dio, 68 - 28887 Omegna (VB) - ITALIA  
Tel. +39.0323.828.004 - Fax +39.0323.828.914  
[www.resista.it](http://www.resista.it) - [info@resista.it](mailto:info@resista.it)

REV 06 - 21/04/2024

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# © calibrated surgical drills