

INTERNAL HEXAGON

NARROW CONNECTION SHORT HEXAGON

NARROW PROSTHETICS



narrow prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato
Healing Abutment - Titanium Gr. 5 - Anodized

IP 3HA3 H 3 mm - Small Gingival Profile
IP 3HA5 H 5 mm - Small Gingival Profile


DM 1/2 Standard Hexagonal Driver



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato
Laboratory Analog - Titanium Gr. 5 - Anodized

IP 3AL H 12 mm - Universal Shape

10Ncm 



Moncone Calcinabile - Acrilico
Castable Abutment - Acrylic

IP 3CA Antirotation Hex - Small Gingival Profile
IP 3CR Rotating - Small Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver

32Ncm 




Moncone Titanio Diritto - Titanio Gr. 5 - PVD TiN
Straight Titanium Abutment - Titanium Gr. 5 - PVD TiN

IP 3ST2 0° - Straight - Small Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver

32Ncm 



Moncone Titanio Angolato - Titanio Gr. 5 - PVD TiN
Angled Titanium Abutment - Titanium Gr. 5 - PVD TiN

IP 3AT15 15° - Angled - Small Gingival Profile
IP 3AT25 25° - Angled - Small Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver

narrow prosthetics

Moncone Provvisorio Estetico - PEEK Mediale

Aesthetic Temporary Abutment - Medical PEEK

IP 3PA Antirotation Hex - Small Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver



Pilastro Temporaneo - Titanio Gr. 5

Temporary Abutment - Titanium Gr. 5

IP 3TA Antirotation Hex - Small Gingival Profile

IP 3TR Rotating - Small Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver



Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Ball Abutment - Titanium Gr. 5 - PVD TiN

IP 3BA1 H 1 mm - Small Gingival Profile - NORMO

IP 3BA2 H 2 mm - Small Gingival Profile - NORMO

IP 3BA4 H 4 mm - Small Gingival Profile - NORMO

DM 1/2 Standard Hexagonal Driver



Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Ball Abutment - Titanium Gr. 5 - PVD TiN

IP 3EQ1 H 1 mm - Small Gingival Profile

IP 3EQ2 H 2 mm - Small Gingival Profile

IP 3EQ3 H 3 mm - Small Gingival Profile

IP 3EQ4 H 4 mm - Small Gingival Profile

IP 3EQ5 H 5 mm - Small Gingival Profile

IP 3EQ6 H 6 mm - Small Gingival Profile



Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

IP X35 H 12mm - Internal Hexagon

DM 1/2 Standard Hexagonal Driver



toronto narrow prosthetics

32Ncm 



Pilastro Toronto 4 Diritto - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

IP 3MS2 H 2 mm - Small Gingival Profile

IP 3MS3 H 3 mm - Small Gingival Profile

DCDM Toronto Straight Abutment Driver

32Ncm 



Pilastro Toronto 4 Angolato - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

IP 3MA173 17° - H 3 mm - Small Gingival Profile

IP 3MA304 30° - H 4 mm - Small Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 Toronto Abutment Torx Driver

15Ncm 



Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

EP MT Rotating - Bonding Technique - Incollaggio Passivo

EP MT-S Rotating - Welding Technique - Saldatura Intraorale

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 Toronto Abutment Torx Driver

10Ncm 



15Ncm 

Pilastri Accessori Toronto - Acrilico

Toronto Accessory Abutments - Acrylic

EP MC Rotating Castable Abutment - Moncone Calcinabile

EP MHA Healing Abutment - Cappetta di Guarigione

DMX 1/2 Toronto Abutment Torx Driver

15Ncm 



Componenti Toronto - Per Impronta e Modello in Gesso

Toronto Components - For Impression and Gypsum Model

EP MI Toronto Impression Abutment - Transfer Impronta

EP MAL Toronto Analog - Analogo da Laboratorio


DMX 1/2 Toronto Abutment Torx Driver

Transfer per Impronta Digitale - Titanio Gr. 5 + Acidificazione

Scan Abutment - Titanium Gr. 5 + Acid Etching

IP 3SA DS Antirotation Hex - **Narrow** Gingival Profile

DM 1/2 Standard Hexagonal Driver



Scansione

Matematiche

Matching delle Matematiche

Link Base per Incollaggio - Titanio Gr. 5

Bonding Technique Base Link - Titanium Gr. 5

IP 3TB DS Antirotation Hex - Adjustable in 3 heights

H1 = 4.0mm - H2 = 5.5mm - H3 = 7.0mm

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver



Link Base per Incollaggio - Titanio Gr. 5

Bonding Technique Base Link - Titanium Gr. 5

IP 3TBR DS Rotating - Adjustable in 3 heights

H1 = 4.0mm - H2 = 5.5mm - H3 = 7.0mm

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver



Analogo per Laboratorio - Titanio Gr. 5

Laboratory Analog - Titanium Gr. 5

IP 3AL DS H 10 mm - Ø 4 Universal Shape

