

EXTERNAL HEXAGON

4.1 REGULAR CONNECTION

UNIVERSAL PROSTHETICS



4.1 universal prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato
Healing Abutment - Titanium Gr. 5 - Anodized

EP 4HA2 H 2 mm - Standard Gingival Profile
EP 4HA4 H 4 mm - Standard Gingival Profile
EP 4HA6 H 6 mm - Standard Gingival Profile

DM 1/2 [Standard Hexagonal Driver](#)

20Ncm 



Pilastro di Guarigione Allargato - Titanio Gr. 5 - Anodizzato
Large Healing Abutment - Titanium Gr. 5 - Anodized

EP 5HA3 H 3 mm - **Large Gingival Profile**
EP 5HA5 H 5 mm - **Large Gingival Profile**


DM 1/2 [Standard Hexagonal Driver](#)



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato
Laboratory Analog - Titanium Gr. 5 - Anodized

EP 4AL H 12 mm - Universal Shape

10Ncm 



Moncone Calcinabile - Acrilico
Castable Abutment - Acrylic

EP 4CA Antirotaion Hex - Standard Gingival Profile
EP 4CR Rotating - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 [Standard Hexagonal Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN
Titanium Abutment - Titanio Gr. 5 - PVD TiN

EP 4ST2 0° - Straight - Standard Gingival Profile
EP 5ST3 0° - Straight - **Large Gingival Profile**

EP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 [Standard Hexagonal Driver](#)

4.1 universal prosthetics

Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanio Gr. 5 - PVD TiN

EP 4AT15 15° - Angled - Standard Gingival Profile

EP 4AT25 25° - Angled - Standard Gingival Profile

EP 5AT15 15° - Angled - **Large Gingival Profile**

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



Moncone UCLA Oro - Base Oro/Platino/Iridio

UCLA Gold Abutment - Gold/Platinum/Iridium Base

EP 4GA Antirotation Hex - Standard Gingival Profile

EP 5GA Antirotation Hex - **Large Gingival Profile**

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



Moncone Provvisorio Estetico - PEEK MedicaLe

Aesthetic Temporary Abutment - Medical PEEK

EP 4PA Antirotation Hex - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 25Ncm



Pilastro Temporaneo Antirotazione - Titanio Gr. 5

Antirotation Temporary Abutment - Titanium Gr. 5

EP 4TA Antirotation Hex - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



Pilastro Temporaneo Rotante - Titanio Gr. 5

Rotating Temporary Abutment - Titanium Gr. 5

EP 4TR Rotating - Standard Gingival Profile

EP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 [Standard Hexagonal Driver](#)



 32Ncm



4.1 universal prosthetics

32Ncm 



Link per Incollaggio TiB - Titanio Gr. 5

Bonding Technique Link TiB - Titanium Gr. 5

EP 4STB Link in Titanio Antirrotazione - Compatibile TiB

EP VIT Vite di Ricambio - Spare Screw
DM 1/2 Standard Hexagonal Driver

32Ncm 



Moncone Estetico - Zirconia

Esthetic Abutment - Zirconium

EP 4SZ 0° - Straight - Standard Gingival Profile

EP 4AZ15 15° - Angled - Standard Gingival Profile

EP VIT Vite di Ricambio - Spare Screw
DM 1/2 Standard Hexagonal Driver

32Ncm 



Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Overdenture Abutment - Titanium Gr. 5 - PVD TiN

EP 4BA1 H 1 mm - Standard Gingival Profile - NORMO

EP 4BA2 H 2 mm - Standard Gingival Profile - NORMO

EP 4BA3 H 3 mm - Standard Gingival Profile - NORMO

EP 4BA4 H 4 mm - Standard Gingival Profile - NORMO

EP 4BA6 H 6 mm - Standard Gingival Profile - NORMO

DM Standard Hexagonal Driver

32Ncm 



Pilastro per Overdenture Equator - Titanio Gr. 5 - PVD TiN

Equator Overdenture Abutment - Titanium Gr. 5 - PVD TiN

EP 4EQ2 H 2 mm - Narrow Gingival Profile


EP 4EQ3 H 3 mm - Narrow Gingival Profile

EP 4EQ4 H 4 mm - Narrow Gingival Profile

EP 4EQ5 H 5 mm - Narrow Gingival Profile

EP 4EQ6 H 6 mm - Narrow Gingival Profile

DMEQ Squared Equator Driver




Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

EP X40 H 12mm - Standard Gingival Profile

EP VITL Vite di Ricambio - Spare Screw
DM 1/2 Standard Hexagonal Driver

toronto 4.1 universal prosthetics



Pilastro Toronto 4 Dritto - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

- EP MS1 H 1 mm - Standard Gingival Profile
- EP MS2 H 2 mm - Standard Gingival Profile
- EP MS3 H 3 mm - Standard Gingival Profile

DCDM Toronto Straight Abutment Driver



Pilastro Toronto 4 Angolato - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

- EP MA172 17° - H 2 mm - Standard Gingival Profile
- EP MA173 17° - H 3 mm - Standard Gingival Profile
- EP MA304 30° - H 4 mm - Standard Gingival Profile

EP VITMA Vite di Ricambio Toronto - Toronto Spare Screw
DMX 1/2 Toronto Abutment Torx Driver



Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

- EP MT Rotating - Bonding Technique - Incollaggio Passivo
- EP MT-S Rotating - Welding Technique - Saldatura Intraorale
- EP MV Vite di Ricambio Toronto - Toronto Spare Screw
- DMX 1/2** Toronto Abutment Torx Driver



toronto 4.1 universal prosthetics

10Ncm



Pilastro Calcinabile Toronto - Acrilico

Toronto Castable Abutment - Acrylic

EP MC Rotating - Standard Gingival Profile

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



15Ncm



Transfer Toronto - Titanio Gr. 5

Toronto Impression Coping - Titanium Gr. 5

EP MI Rotating - Standard Gingival Profile

EP MVL Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)



Analogo Toronto - Titanio Gr. 5

Toronto Analog - Titanium Gr. 5

EP MAL H 10 mm - Universal Shape

15Ncm



Pilastro di Guarigione Toronto - PE

Toronto Healing Abutment - PE

EP MHA H 4 mm - Standard Gingival Profile

EP MHA6 H 6 mm - Standard Gingival Profile

DMX 1/2 [Toronto Abutment Torx Driver](#)

15Ncm



15Ncm



Viti di Ricambio Toronto - Titanio Gr. 5

Toronto Spare Screws - Titanium Gr. 5

EP MV Prosthetic Spare Screw

EP MVL Impression Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)