

narrow platform prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato
Healing Abutment - Titanium Gr. 5 - Anodized

SP 3HA3 H 3 mm - Narrow Gingival Profile
SP 3HA5 H 5 mm - Narrow Gingival Profile


DMX 1/2 [Toronto Abutment Torx Driver](#)



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato
Laboratory Analog - Titanium Gr. 5 - Anodized

SP 3AL H 12 mm - Narrow Shape

10Ncm 



Moncone Calcinabile - Acrilico
Castable Abutment - Acrylic

SP 3CA Antirotaion Hex - Narrow Gingival Profile
SP 3CR Rotating - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw
DMX 1/2 [Toronto Abutment Torx Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN
Titanium Abutment - Titanio Gr. 5 - PVD TiN

SP 3ST 0° - Straight - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw
DMX 1/2 [Toronto Abutment Torx Driver](#)

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN
Titanium Abutment - Titanio Gr. 5 - PVD TiN

SP 3AT15 15° - Angled - Narrow Gingival Profile
SP 3AT25 25° - Angled - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw
DMX 1/2 [Toronto Abutment Torx Driver](#)

narrow platform prosthetics

Pilastro Temporaneo - Titanio Gr. 5

Temporary Abutment - Titanium Gr. 5

SP 3TA Antirotaion Hex - Narrow Gingival Profile

SP 3TR Rotating - Narrow Gingival Profile

SP VIT33 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

SP 3IA H 12mm - Narrow Gingival Profile

SP VITL33 Vite di Ricambio Lunga - Long Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



regular platform prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

SP 4HA3 H 3 mm - Standard Gingival Profile

SP 4HA5 H 5 mm - Standard Gingival Profile

DMX 1/2

Toronto Abutment Torx Driver



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

SP 4AL H 12 mm - Standard Shape

10Ncm 



Moncone Calcinabile Antirotazione - Acrilico

Antirootation Castable Abutment - Acrylic

SP 4CA Antirootation Hex - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

10Ncm 



Moncone Calcinabile Rotante - Acrilico

Rotating Castable Abutment - Acrylic

SP 4CR Rotating - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

32Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanio Gr. 5 - PVD TiN

SP 4ST 0° - Straight - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

regular platform prosthetics

Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanio Gr. 5 - PVD TiN

SP 4AT15 15° - Angled - Standard Gingival Profile

SP 4AT25 25° - Angled - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Link per Incollaggio TiB - Titanio Gr. 5

Bonding Technique Link TiB - Titanium Gr. 5

EP 4STB Link in Titanio Antirotazione - Compatibile TiB

EP VIT Vite di Ricambio - Spare Screw

DM 1/2 Standard Hexagonal Driver



Pilastro Temporaneo Antirotazione - Titanio Gr. 5

Antirotation Temporary Abutment - Titanium Gr. 5

SP 4TA Antirotation Hex - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Pilastro Temporaneo Rotante - Titanio Gr. 5

Rotating Temporary Abutment - Titanium Gr. 5

SP 4TR Rotating - Standard Gingival Profile

SP VIT40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver




Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5


SP 4IA H 12mm - Standard Gingival Profile

SP VITL40 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



toronto universal prosthetics



Pilastro Toronto 4 Dritto - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

- EP MS2 H 2 mm - Standard Gingival Profile
- EP MS3 H 3 mm - Standard Gingival Profile

DCDM Toronto Straight Abutment Driver




Pilastro Toronto 4 Angolato - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

- EP MA172 17° - H 2 mm - Standard Gingival Profile
- EP MA173 17° - H 3 mm - Standard Gingival Profile
- EP MA304 30° - H 4 mm - Standard Gingival Profile

EP VITMA Vite di Ricambio Toronto - Toronto Spare Screw
DMX 1/2 Toronto Abutment Torx Driver



Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

- EP MT Rotating - Bonding Technique - Incollaggio Passivo
- EP MT-S Rotating - Welding Technique - Saldatura Intraorale
- EP MV Vite di Ricambio Toronto - Toronto Spare Screw
- DMX 1/2** Toronto Abutment Torx Driver



toronto universal prosthetics

10Ncm



Pilastro Calcinabile Toronto - Acrilico

Toronto Castable Abutment - Acrylic

EP MC Rotating - Standard Gingival Profile

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)

15Ncm



Pilastro di Guarigione Toronto - PE

Toronto Healing Abutment - PE

EP MHA H 4 mm - Standard Gingival Profile

EP MHA6 H 6 mm - Standard Gingival Profile

DMX 1/2 [Toronto Abutment Torx Driver](#)



15Ncm



Transfer Toronto - Titanio Gr. 5

Toronto Impression Coping - Titanium Gr. 5

EP MI Rotating - Standard Gingival Profile

DMX 1/2 [Toronto Abutment Torx Driver](#)



Analogo Toronto - Titanio Gr. 5

Toronto Analog - Titanium Gr. 5

EP MAL H 10 mm - Universal Shape

15Ncm



15Ncm



Viti di Ricambio Toronto - Titanio Gr. 5

Toronto Spare Screws - Titanium Gr. 5

EP MV Prosthetic Spare Screw

EP MVL Impression Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)

wide platform prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato

Healing Abutment - Titanium Gr. 5 - Anodized

SP 5HA3 H 3 mm - Large Gingival Profile

SP 5HA5 H 5 mm - Large Gingival Profile

DMX 1/2

Toronto Abutment Torx Driver

Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato

Laboratory Analog - Titanium Gr. 5 - Anodized

SP 5AL H 12 mm - Large Shape



Moncone Calcinabile - Acrilico

Castable Abutment - Acrylic

SP 5CA Antirotação Hex - Narrow Gingival Profile

SP 5CR Rotating - Large Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

10Ncm 



Moncone Titanio - Titanio Gr. 5 - PVD TiN

Titanium Abutment - Titanium Gr. 5 - PVD TiN

SP 5ST 0° - Straight - Large Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

32Ncm 



Pilastro Temporaneo con Esagono - Titanio Gr. 5

Antirotação Hex Temporary Abutment - Titanium Gr. 5

SP 5TA Antirotação Hex - Narrow Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2

Toronto Abutment Torx Driver

32Ncm 



wide platform prosthetics

Pilastro Temporaneo senza Esagono - Titanio Gr. 5

Rotating Temporary Abutment - Titanium Gr. 5

SP 5TR Rotating - Narrow Gingival Profile

SP VIT50 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver



Transfer per Impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

SP 5IA H 12mm - Narrow Gingival Profile

SP VITL50 Vite di Ricambio - Spare Screw

DMX 1/2 Toronto Abutment Torx Driver

