

# INTERNAL HEXAGON NARROW CONNECTION

## NARROW PROSTHETICS



# narrow prosthetics

20Ncm 



**Pilastro di Guarigione** - Titanio Gr. 5 - Anodizzato  
**Healing Abutment** - Titanium Gr. 5 - Anodized

**IP 3HA3**      H 3 mm - Small Gingival Profile  
**IP 3HA5**      H 5 mm - Small Gingival Profile


**DM 1/2**      [Standard Hexagonal Driver](#)



**Analogo per Laboratorio** - Titanio Gr. 5 - Anodizzato  
**Laboratory Analog** - Titanium Gr. 5 - Anodized

**IP 3AL**      H 12 mm - Universal Shape

10Ncm 



**Moncone Calcinabile** - Acrilico  
**Castable Abutment** - Acrylic

**IP 3CA**      Antirotation Hex - Small Gingival Profile  
**IP 3CR**      Rotating - Small Gingival Profile

**IP VITG**      Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2**      [Standard Hexagonal Driver](#)

32Ncm 




**Moncone Titanio Diritto** - Titanio Gr. 5 - PVD TiN  
**Straight Titanium Abutment** - Titanium Gr. 5 - PVD TiN

**IP 3ST2**      0° - Straight - Small Gingival Profile

**IP VITG**      Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2**      [Standard Hexagonal Driver](#)

32Ncm 



**Moncone Titanio Angolato** - Titanio Gr. 5 - PVD TiN  
**Angled Titanium Abutment** - Titanium Gr. 5 - PVD TiN

**IP 3AT15**      15° - Angled - Small Gingival Profile  
**IP 3AT25**      25° - Angled - Small Gingival Profile

**IP VITG**      Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2**      [Standard Hexagonal Driver](#)

# narrow prosthetics

## Moncone Provvisorio Estetico - PEEK Mediale

Aesthetic Temporary Abutment - Medical PEEK

**IP 3PA** Antirotation Hex - Small Gingival Profile

**IP VITG** Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2** Standard Hexagonal Driver



## Pilastro Temporaneo - Titanio Gr. 5

Temporary Abutment - Titanium Gr. 5

**IP 3TA** Antirotation Hex - Small Gingival Profile

**IP 3TR** Rotating - Small Gingival Profile

**IP VITG** Vite di Ricambio Gold - Gold Spare Screw  
**DM 1/2** Standard Hexagonal Driver



## Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Ball Abutment - Titanium Gr. 5 - PVD TiN

**IP 3BA1** H 1 mm - Small Gingival Profile - NORMO

**IP 3BA2** H 2 mm - Small Gingival Profile - NORMO

**IP 3BA4** H 4 mm - Small Gingival Profile - NORMO

**DM 1/2** Standard Hexagonal Driver



## Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Ball Abutment - Titanium Gr. 5 - PVD TiN

**IP 3EQ1** H 1 mm - Small Gingival Profile

**IP 3EQ2** H 2 mm - Small Gingival Profile

**IP 3EQ3** H 3 mm - Small Gingival Profile

**IP 3EQ4** H 4 mm - Small Gingival Profile

**IP 3EQ5** H 5 mm - Small Gingival Profile

**IP 3EQ6** H 6 mm - Small Gingival Profile



## Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

**IP X35** H 12mm - Internal Hexagon

**DM 1/2** Standard Hexagonal Driver



# toronto narrow prosthetics

32Ncm 



**Pilastro Toronto 4 Diritto** - Titanio Gr. 5

Toronto 4 Straight Abutment - Titanium Gr. 5

IP 3MS2 H 2 mm - Small Gingival Profile

IP 3MS3 H 3 mm - Small Gingival Profile

**DCDM** Toronto Straight Abutment Driver

32Ncm 



**Pilastro Toronto 4 Angolato** - Titanio Gr. 5

Toronto 4 Angled Abutment - Titanium Gr. 5

IP 3MA173 17° - H 3 mm - Small Gingival Profile

IP 3MA304 30° - H 4 mm - Small Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver

15Ncm 



**Pilastro Temporaneo Toronto** - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5

EP MT Rotating - Bonding Technique - Incollaggio Passivo

EP MT-S Rotating - Welding Technique - Saldatura Intraorale

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

**DMX 1/2** Toronto Abutment Torx Driver

10Ncm 



15Ncm 



**Pilastri Accessori Toronto** - Acrilico

Toronto Accessory Abutments - Acrylic

EP MC Rotating Castable Abutment - Moncone Calcinabile

EP MHA Healing Abutment - Cappetta di Guarigione

**DMX 1/2** Toronto Abutment Torx Driver

15Ncm 



**Componenti Toronto** - Per Impronta e Modello in Gesso

Toronto Components - For Impression and Gypsum Model

EP MI Toronto Impression Abutment - Transfer Impronta

EP MAL Toronto Analog - Analogo da Laboratorio


**DMX 1/2** Toronto Abutment Torx Driver

**Transfer per Impronta Digitale - Titanio Gr. 5 + Acidificazione**

Scan Abutment - Titanium Gr. 5 + Acid Etching

IP 3SA DS      Antirotation Hex - **Narrow** Gingival Profile

DM 1/2      Standard Hexagonal Driver



**Link Base per Incollaggio - Titanio Gr. 5**

Bonding Technique Base Link - Titanium Gr. 5

IP 3TB DS      Antirotation Hex - Adjustable in 3 heights

H1 = 4.0mm - H2 = 5.5mm - H3 = 7.0mm

IP VITG      Vite di Ricambio Gold - Gold Spare Screw

DM 1/2      Standard Hexagonal Driver



**Link Base per Incollaggio - Titanio Gr. 5**


Bonding Technique Base Link - Titanium Gr. 5

IP 3TBR DS      Rotating - Adjustable in 3 heights

H1 = 4.0mm - H2 = 5.5mm - H3 = 7.0mm

IP VITG      Vite di Ricambio Gold - Gold Spare Screw

DM 1/2      Standard Hexagonal Driver



**Analogo per Laboratorio - Titanio Gr. 5**

Laboratory Analog - Titanium Gr. 5

IP 3AL DS      H 10 mm - Ø 4 Universal Shape

