

wide prosthetics

20Ncm 



Pilastro di Guarigione - Titanio Gr. 5 - Anodizzato
Healing Abutment - Titanium Gr. 5 - Anodized

IP 5HA3 H 3 mm - Large Gingival Profile

IP 5HA5 H 5 mm - Large Gingival Profile


DM 1/2 Standard Hexagonal Driver



Analogo per Laboratorio - Titanio Gr. 5 - Anodizzato
Laboratory Analog - Titanium Gr. 5 - Anodized

IP 5AL H 12 mm - Universal Shape

10Ncm 



Moncone Calcinabile - Acrilico
Castable Abutment - Acrylic


IP 5CA Antirotation Hex - Large Gingival Profile

IP 5CR Rotating - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver

32Ncm 




Moncone Titanio Diritto - Titanio Gr. 5 - PVD TiN
Straight Titanium Abutment - Titanium Gr. 5 - PVD TiN

IP 5ST3 0° - Straight - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver

32Ncm 



Moncone Titanio Angolato - Titanio Gr. 5 - PVD TiN
Angled Titanium Abutment - Titanium Gr. 5 - PVD TiN

IP 5AT15 15° - Angled - Large Gingival Profile

IP 5AT25 25° - Angled - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver

wide prosthetics



Pilastro Temporaneo da Incollaggio - Titanio Gr. 5

Temporary Glueing Abutment - Titanium Gr. 5

IP 4TA-W Antirotation Hex - **Waisted** Gingival Profile CPS

IP 4TR-W Rotating - **Waisted** Gingival Profile CPS

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver




Pilastro Temporaneo da Saldatura - Titanio Gr. 5

Temporary Welding Abutment - Titanium Gr. 5

IP 4TR-S Rotating - **Waisted** Gingival Profile CPS

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver




Moncone Ucla Oro - Base Oro/Platino/Iridio

Ucla Gold Abutment - Gold/Platinum/Iridium Base

IP 5GA Antirotation Hex - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver



Moncone UCLA Cobalto Cromo - CoCr


UCLA Abutment Cobalt Chrome - CoCr

IP 4CCA Antirotation Hex - Standard Gingival Profile

IP 4CCA-W Antirotation Hex - **Waisted** Gingival Profile CPS

IP 4CCR Rotating Hex - Standard Gingival Profile

IP 4CCR-W Rotating Hex - **Waisted** Gingival Profile CPS




Moncone Provvisorio Estetico - PEEK Medica

Aesthetic Temporary Abutment - Medical PEEK

IP 5PA Antirotation Hex - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw

DM 1/2 Standard Hexagonal Driver



wide prosthetics

32Ncm 




Moncone Estetico Diritto - Zirconia + Link in Titanio Gr. 5

Straight Titanium Abutment - Zirconium + Titanium Gr. 5 Link

IP 5SZ 0° - Straight - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver

32Ncm 



Moncone Estetico Angolato - Zirconia + Link in Titanio Gr. 5

Angled Titanium Abutment - Zirconium + Titanium Gr. 5 Link

IP 5AZ15 15° - Angled - Large Gingival Profile

IP VITG Vite di Ricambio Gold - Gold Spare Screw
DM 1/2 Standard Hexagonal Driver

32Ncm 



Pilastro per Overdenture - Titanio Gr. 5 - PVD TiN

Overdenture Abutment - Titanium Gr. 5 - PVD TiN

IP 4BA1-W H 1 mm - **Waisted** Gingival Profile CPS - NORMO

IP 4BA2-W H 2 mm - **Waisted** Gingival Profile CPS - NORMO

IP 4BA3-W H 3 mm - **Waisted** Gingival Profile CPS - NORMO

IP 4BA4-W H 4 mm - **Waisted** Gingival Profile CPS - NORMO

IP 4BA6-W H 6 mm - **Waisted** Gingival Profile CPS - NORMO

DM Standard Hexagonal Driver

32Ncm 



Pilastro per Overdenture Equator - Titanio Gr. 5 - PVD TiN

Equator Overdenture Abutment - Titanium Gr. 5 - PVD TiN

IP 4EQ1 H 1 mm - Narrow Gingival Profile

IP 4EQ2 H 2 mm - Narrow Gingival Profile

IP 4EQ3 H 3 mm - Narrow Gingival Profile

IP 4EQ4 H 4 mm - Narrow Gingival Profile

IP 4EQ5 H 5 mm - Narrow Gingival Profile

DMEQ Squared Equator Driver



 25Ncm

Transfer per impronta Pick Up - Titanio Gr. 5

Pick Up Impression Coping - Titanium Gr. 5

IP X50 H 12mm - Internal Hexagon

DM 1/2 Standard Hexagonal Driver

Pilastro Toronto Dritto - Titanio Gr. 5

Toronto Straight Abutment - Titanium Gr. 5

IP MS1-W H 1 mm - Standard Gingival Profile

IP MS2-W H 2 mm - Standard Gingival Profile

IP MS3-W H 3 mm - Standard Gingival Profile

DCDM [Toronto Straight Abutment Driver](#)

 32Ncm



Pilastro Toronto Angolato 17° - Titanio Gr. 5

17° Toronto Angled Abutment - Titanium Gr. 5

IP MA172-W 17° - H 2 mm - Standard Gingival Profile

IP MA173-W 17° - H 3 mm - Standard Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)

 32Ncm



Pilastro Toronto Angolato 30° - Titanio Gr. 5

30° Toronto Angled Abutment - Titanium Gr. 5

IP MA303-W 30° - H 3 mm - Standard Gingival Profile

IP MA304-W 30° - H 4 mm - Standard Gingival Profile

IP MA305-W 30° - H 5 mm - Standard Gingival Profile

IP VITMA Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)

 32Ncm



Pilastro Temporaneo Toronto - Titanio Gr. 5

Toronto Temporary Abutment - Titanium Gr. 5


EP MT Rotating - Bonding Technique - Incollaggio Passivo

EP MT-S Rotating - Welding Technique - Saldatura Intraorale

EP MV Vite di Ricambio Toronto - Toronto Spare Screw

DMX 1/2 [Toronto Abutment Torx Driver](#)

15Ncm



Pilastri di Gyuarigione e Transfer Toronto

Toronto Healing Abutment and Impression Coping

EP MI Rotating - Standard Gingival Profile

EP MHA H 4 mm - Standard Gingival Profile

EP MHA6 H 6 mm - Standard Gingival Profile

DMX 1/2 [Toronto Abutment Torx Driver](#)

15Ncm

 15Ncm

